## EurobioPlex

### *SARS-CoV-2 SNPs*

## RT-PCR EN TEMPS REEL

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Reference EBX-046 Version 5.00 – July 2021

### **Validated on:**

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)

Instructions for use

Available on [www.eurobio-scientific.com](http://www.eurobio-scientific.com)

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## Table of contents

|                                                                                                              |    |
|--------------------------------------------------------------------------------------------------------------|----|
| 1. Introduction.....                                                                                         | 30 |
| 2. Purpose of the system.....                                                                                | 31 |
| 3. Symbols .....                                                                                             | 32 |
| 4. Principle of detection .....                                                                              | 33 |
| 5. Content of the kit .....                                                                                  | 34 |
| 6. Storage.....                                                                                              | 35 |
| 7. Materials required not provided .....                                                                     | 35 |
| 8. Real-time PCR instrument .....                                                                            | 35 |
| 9. Cautions and note .....                                                                                   | 36 |
| 10. Procedure .....                                                                                          | 37 |
| 11. Validation of the experiment .....                                                                       | 41 |
| 12. Data analysis and interpretation .....                                                                   | 44 |
| 13. Performance analysis .....                                                                               | 45 |
| 14. Comparative analysis of the EBX-046 kit's diagnostic performance versus a commercial reference kit ..... | 52 |
| 15. Bibliography.....                                                                                        | 53 |
| 16. Quality control.....                                                                                     | 53 |
| 17. Waste disposal.....                                                                                      | 54 |
| 18. Incident report.....                                                                                     | 54 |
| 19. Technical assistance .....                                                                               | 54 |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## 1. Introduction

The SARS-CoV-2 virus emerged at the end of 2019 in the city of Wuhan in China. It is classified as a member of the Coronaviridae family and the Sarbecovirus subgroup. Its genome consists of 29903 ribonucleic acid (RNA) bases. SARS-CoV-2 is the 7th identified coronavirus infecting humans, after human coronaviruses (HCoV), SARS-CoV (coronavirus causing severe acute respiratory syndrome) and MERS-CoV (coronavirus inducing Middle East respiratory syndrome).

The sequences of SARS-CoV-2 show similarities to those of the betacoronaviruses found in bats. SARS-CoV-2 is genetically distinct from other human coronaviruses.

The clinical profile is variable, ranging from flu-like symptoms, fever, cough, difficulty breathing, to pneumonia, to a severe respiratory syndrome that can be fatal. The lethality rate of the virus is estimated to be between 2-3%. SARS-CoV-2 is highly contagious with an estimated 173 million confirmed cases worldwide since early June 2021.

Several variants and mutations of SARS-CoV-2 have demonstrated an impact on public health (increased transmissibility, severity of infection or immune escape). They emerged at the end of 2020 in distinct geographical areas and are designated as variants of concern (VOC):

- Alpha variant from the United Kingdom (UK): lineage B.1.1.7 with mutation N501Y.
- Beta variant from South Africa (SA): lineage B.1.351 with mutation E484K.
- Gamma variant from Brazil (BR): lineage P.1 with mutation E484K.
- Delta variant from India: lineage B.1.617.2 with mutation L452R.
- Kappa variant from India: lineage B.1.617.1 with the mutations L452R and E484Q.
- The L452R mutation is also found in other lineages, including A.27, the origin of which is unclear, and C.36, detected in Egypt. An L452Q variation of lineage C.37 has also been detected in Peru.

In May 2021, 80% of Covid-19 infections were reported with the UK B.1.1.7 variant. This variant is recognised and the current vaccines on the market provide protection against it. The other registered cases were related to the other circulating variants. The interest in the diagnosis of the distinct Covid-19 variants is therefore to focus on the recent E484K, E484Q and L452R mutations for which currently available vaccines may have limited or no effect, being mutations responsible for immune escape. The new variant identification strategy aims to identify immune escape mutations as a priority, rather than distinguishing variants. This strategy will enable patients with variants with this type of mutation to be treated more quickly.

The EurobioPlex SARS-CoV-2 SNPs kit was designed by aligning in silico all the SARS-CoV-2 sequences known to date.

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## 2. Purpose of the system

The EurobioPlex SARS-CoV-2 SNPs is a real-time polymerase chain reaction (RT-PCR) assay that, in a single test, identifies the presence of SARS-CoV-2 and distinguishes the E484K, E484Q and L452R mutations in a qualitative manner.

The RNA extract is the starting material for the EurobioPlex SARS-CoV-2 kit. It is up to the user to use extraction methods adapted to the tested samples. The kit has been validated on nasopharyngeal samples.

The EurobioPlex SARS-CoV-2 SNPs kit is an *in Vitro Diagnostic Medical Device* for the detection of SARS-CoV-2 and mutations of interest. This kit is available for first-line diagnostic purposes.

The EurobioPlex SARS-CoV-2 SNPs test must be used by qualified medical laboratory staff. It is for single use only and should not be reused after usage.

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Manufacturer

CE marked product

In Vitro Diagnostic Medical Device

Instructions for use

Single use

Caution

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

## 4. Principle of detection

The EurobioPlex SARS-CoV-2 SNPs kit is a qualitative ribonucleic acid (RNA) amplification assay for the analysis of:

- The presence of the two targets (RdRp1/RdRp2) of SARS-CoV-2 on the RdRp gene, in the same channel.
- The presence of E484K, E484Q and L452R mutations and their distinction in three distinct detection channels.
- The presence of a human housekeeping gene used as a control for the quality of the sample, for RNA extraction and RT-PCR inhibition. The control is used to identify potential variations that may occur during the RNA extraction from biological samples and real-time RT-PCR amplification steps. This ensures that a negative result cannot be due to incorrect RNA extraction and/or the presence of an excessive quantity of RT-PCR inhibitors.

The RNA extract is the source material for the EurobioPlex SARS-CoV-2 SNPs kit. The RNA is extracted from a patient's nasopharyngeal sample.

SARS-CoV-2 RNA is detected using specific probes for each mutation or gene labelled as FAM, HEX, Texas Red, Quasar 705. The endogenous internal control is detected using a Cy5 labelled probe (Table 1). During elongation of the amplification product, the probes emit specific fluorescence following their hydrolysis. The measurement of the fluorescence intensity in real time, emitted and measured individually by an optical system during the PCR, is related to the accumulation of the amplification products.

**Table 1: Detection of targets by fluorophores**

| Targets                     | Fluorophore | Excitation | Emission |
|-----------------------------|-------------|------------|----------|
| 484K                        | FAM         | 495 nm     | 515 nm   |
| 484Q                        | HEX         | 535 nm     | 555 nm   |
| RdRp1/RdRp2                 | Texas Red   | 585 nm     | 605 nm   |
| Endogenous internal control | Cy5         | 650 nm     | 670 nm   |
| 452R                        | Quasar 705  | 683 nm     | 703 nm   |

Recommended equivalent channels on different PCR instruments:

- **FAM channel** (CFX96, QuantStudio 5, QuantStudio 6), Green channel (Rotor-Gene Q),
- **HEX channel** (CFX96), VIC channel (QuantStudio 5, QuantStudio 6), Yellow channel (Rotor-Gene Q),

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

- **Texas Red channel** (CFX96, QuantStudio 6), Orange channel (Rotor-Gene Q), Rox channel (QuantStudio 5),
- **Cy5 channel** (CFX96, QuantStudio 5, QuantStudio 6), Red channel (Rotor-Gene Q),
- **Quasar 705 channel** (CFX96, QuantStudio 5, QuantStudio 6), Purple channel (Rotor-Gene Q).

## 5. Content of the kit

The Eurobioplex SARS-CoV-2 SNPs real-time RT-PCR kit is ready to use and contains the reagents and enzymes necessary for the detection of the virus (Table 2).

Fluorescence is emitted and measured individually by an optical system during the PCR process. Detection of the amplified fragments is performed by a fluorometer using the channels listed in Table 1.

**Table 2: Content of the kit**

| Cap color          | Content of the kit                   | 48 reactions | 96 reactions | 192 reactions | Reconstitution |
|--------------------|--------------------------------------|--------------|--------------|---------------|----------------|
| <b>Red</b>         | Enzyme Taq polymere-<br>rase         | 650 µL       | 1500 µL      | 2 X 1250 µL   | Ready to use   |
| <b>Brown</b>       | Enzyme Reverse<br>transcriptase      | 15 µL        | 30 µL        | 40 µL         | Ready to use   |
| <b>Transparent</b> | Oligomix                             | 260 µL       | 600 µL       | 1000 µL       | Ready to use   |
| <b>Yellow</b>      | Positive Control (CP)                | 80 µL        | 160 µL       | 320 µL        | Ready to use   |
| <b>Blue</b>        | Water = Negative<br>Control (CN-H2O) | 1 mL         | 1 mL         | 1 mL          | Ready to use   |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## 6. Storage

All reagents should be stored between -15°C and -22°C.

All reagents can be used until the expiration date indicated on the kit label.

**The kit can be used after the initial opening for a maximum of 5 freezing and thawing cycles. Beyond that, the performance of the kit is not guaranteed.**

## 7. Materials required not provided

- ◇ Biological cabinet
- ◇ Real time PCR instrument
- ◇ Centrifuge for microtubes
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR reaction
- ◇ Micropipettes
- ◇ DNase-free and RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (talc-free)

## 8. Real-time PCR instrument

The EurobioPlex SARS-CoV-2 SNPs kit has been developed and validated for use with the following real-time PCR instruments:

- CFX96™ Real Time PCR detection system (Bio-Rad) with analysis on CFX Manager version 3.1 (Bio-Rad)

## 9. Cautions and note

**Read these instructions carefully before starting the procedure.**

- ◇ This experimentation should be performed by medical laboratory technicians.
- ◇ The local and national biosafety regulations in place for the detection of SARS-CoV-2 must be followed strictly at all times, especially in laboratories and with laboratory equipment in agreement.
- ◇ Ensure that instruments have been installed, calibrated and maintained in accordance with the manufacturer's recommendations.
- ◇ It is the responsibility of the user to use equipment other than that validated, and in such cases performance is not guaranteed.
- ◇ Clinical samples should be considered as potentially infectious material and should be prepared in a laminar flow cabinet.
- ◇ This experiment should be performed according to good laboratory practice.
- ◇ Do not use this kit after the stated expiration date.
- ◇ The kit is shipped under dry ice, and the kit components should arrive frozen. If one or more of the components arrive thawed, or if the tubes have been damaged in transit, contact Eurobio Scientific.
- ◇ Avoid freeze/thaw cycles of reagents as this may lead to a decrease in assay sensitivity.
- ◇ Once the reagents have been thawed, centrifuge the tubes briefly before use.
- ◇ The use of ice or an ice pack is recommended in case of long delays due to, for example, a large number of samples to be processed or high temperatures.
- ◇ It is recommended to define three separate work areas: 1) RNA isolation, 2) Preparation of the reaction mixture, and 3) Amplification/Detection of the amplified products.
- ◇ It is recommended that the positive control be opened and handled separately from the biological samples to be tested and the other kit components to avoid cross-contamination.
- ◇ Distinct lab coat and gloves (talco-free) should be worn in each work area.
- ◇ Pipettes, reagents and other working materials should not be moved between these areas.
- ◇ Special care should be taken to maintain the purity of reagents and reaction mixtures.
- ◇ Appropriate RNA preparation/extraction methods for quality RNA production and RT-PCR application should be used, especially to avoid contamination with RNAs.
- ◇ Use filter tips for micropipettes, RNase-free and DNase-free.
- ◇ Do not pipette with the mouth. Do not eat, drink or smoke in the laboratory.
- ◇ Avoid aérosols.

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

## 10. Procedure

### a. Samples collection

- ◇ Collect samples in sterile tubes.
- ◇ It is the responsibility of the user to control their own sample collection, transport, storage and extraction conditions to ensure that RNA extraction using appropriate systems produces RNA of high quality.
- ! ◇ It is recommended that samples are extracted immediately, or stored according to the sample storage recommendations before extraction (Table 3).

**Table 3: Storage recommendations before use**

| Recommendations for maximum storage of samples before extraction |                   |
|------------------------------------------------------------------|-------------------|
| Ambient temperature                                              | 2 h               |
| 4°C                                                              | 72 h              |
| -20°C (preferably -80°C)                                         | Long-term storage |

| Caution                                                                             |                                                                                                                                                                                                                                                                                                                                                                                           |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <ul style="list-style-type: none"><li>◇ The user may refer to the recommendations issued by the World Health Organization or the French National Authority for Health for the proper storage of samples.</li><li>◇ Extracted RNA should be stored at -80°C.</li><li>◇ The transport of clinical samples is subject to local regulations for the transport of infectious agents.</li></ul> |

### b. RNA extraction

It is the responsibility of the user to ensure that the nucleic acid extraction system used is compatible with real-time RT-PCR technology. For this kit, we recommend using virus RNA extraction methods suitable for respiratory specimens, and to refer to the instructions of the supplier of the extraction kit used.

In the EBX-046 kit, the endogenous internal control on the Cy5 channel is already present in the clinical sample. Properly detected after amplification, the endogenous internal control is used to ensure the quality of the sample and the extraction. The performance described in this document was obtained with the following extraction technique: Maelstrom 9600 (TanBEAD) and the OptiPure Viral Auto Plate kit (ref: W665A46).

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## c. Real-time RT-PCR

### General note:

- The positive control contains high concentrations of matrix. Handling must be done carefully to avoid contamination.
- To confirm that the RT-PCR is working correctly, it is necessary to test the positive control, as well as the negative control (water supplied = CN-H<sub>2</sub>O) (see II-2/6 of the real-time RT-PCR protocol).

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

Diagram of the procedure:

## 1 - PREPARATION OF MASTERMIX

| Number of reactions          | N+3*            |
|------------------------------|-----------------|
| Enzyme Taq polymerase        | (N+3) x 12,5 µL |
| Enzyme Reverse transcriptase | (N+3) x 0,2 µL  |
| Water                        | (N+3) x 2,3 µL  |
| Oligomix                     | (N+3) x 5 µL    |
| Total volume Mastermix       | (N+3) x 20 µL   |

## 2 - PREPARATION OF REACTIONS

### Sample

20 µL Mastermix  
+  
5 µL RNA sample

### Positive Control

20 µL Mastermix  
+  
5µL CP

### Negative Control

20 µL Mastermix  
+  
5 µL Molecular biology water (CN-H2O)

## 3 – REAL TIME PCR INSTRUMENT

| Program               | Temperature | Duration | Cycle(s) |                             |
|-----------------------|-------------|----------|----------|-----------------------------|
| Reverse Transcription | 45°C        | 10 min   | 1        | -                           |
| Denaturation          | 95°C        | 3 min    | 1        | -                           |
| Amplification         | 95°C        | 15 sec   | 45       | -                           |
|                       | 58°C        | 30 sec   |          | Acquisition of fluorescence |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

## d. Detailed procedure

- 1) Ensure that the reagents are completely thawed. Homogenise the enzyme tubes, and vortex the Oligomix and CP for about 15 seconds, then centrifuge briefly.
- 2) Prepare the Mastermix as below. N is the number of reactions (including positive and negative controls. Plan to prepare enough reagents for at least N+3 reactions.

| Nombre de réactions          | N+3*            |
|------------------------------|-----------------|
| Enzyme Taq polymérase        | (N+3) x 12,5 µL |
| Enzyme Reverse transcriptase | (N+3) x 0,2 µL  |
| Eau                          | (N+3) x 2,3 µL  |
| Oligomix                     | (N+3) x 5 µL    |
| Volume total Mastermix       | (N+3) x 20 µL   |

*\*For smaller series (≤10): preparing for N+2 is sufficient.*

- 3) Homogenise the Mastermix prepared in 2) and centrifuge briefly
- 4) Dispense 20 µL of Mastermix into separate microplate tubes/wells using a micropipette and filter tips.
- 5) Add 5 µL of extracted RNA sample.
- 6) In parallel perform the following controls:
  - Positive Control:
    - 20 µL of Mastermix + 5 µL of CP.
  - Negative Control:
    - 20 µL of Mastermix + 5 µL supplied water (CN-H<sub>2</sub>O).
- 7) Immediately close the tubes, or plate with an adhesive film to avoid contamination.
- 8) Centrifuge briefly to collect the reaction mixture at the bottom of the tubes or microplate wells.
- 9) Run the following program on the real-time PCR instrument.

| Program               | Temperature | Duration | Cycle(s) |                             |
|-----------------------|-------------|----------|----------|-----------------------------|
| Reverse Transcription | 45°C        | 10 min   | 1        | -                           |
| Denaturation          | 95°C        | 3 min    | 1        | -                           |
| Amplification         | 95°C        | 15 sec   | 45       | -                           |
|                       | 58°C        | 30 sec   |          | Acquisition of fluorescence |

**Note 1:** On the CFX96™ (Bio-Rad), start the run from version 1.6 or later, then analyse with version 3.1 (see § Validation of the experiment). **Please ensure that white PCR plates are used for appropriate reading in all channels.**

**Note 2:** For Applied Biosystems instruments, select "NONE" in "PASSIVE REFERENCE".

**Note 3:** On QuantStudio 5® (Applied Biosystems), start the run from QuantStudio Design & Analysis software v1.5, then analyse with Design & Analysis software v1.5 or v2.5

# Fiche technique EurobioPlex SARS-CoV-2 SNPs


---

**Note 4:** On QuantStudio 6® (Applied Biosystems), start the run from QuantStudio Real-Time PCR software v1.3, then analyze with Design & Analysis software v1.3 or v2.5

## 11. Validation of the experiment


Post-acquisition data analysis on a CFX96 PCR instrument (Bio-Rad) must be performed using version 3.1 of the CFX Manager software (Bio-Rad). In order to change to this version from a run started on an earlier version, please follow the procedure below: At the end of the run, the data file with the extension .pcrd must be opened and processed with CFX Manager (Bio-Rad) version 3.1.

If the run was started with CFX Manager v1.6 for example, to open a data file with CFX Manager v3.1, click on the CFX Manager v3.1 icon. The home screen appears.


Click on "File" and select "Open" then "Data File".

# Fiche technique EurobioPlex SARS-CoV-2 SNPs


Select the file you need to analyse and click on "Open".

The "Drift correction" option must be applied in the "Settings" tab as shown in the image below: click on the "Settings" tab, then on "Baseline Setting" and on "Apply Fluorescence Drift Correction".


To optimise the run analysis, select the "Log Scale" checkbox for each channel analysed. Then place the threshold bar above the background noise corresponding to the middle of the exponential phase. The 5 channels of interest have high RFUs and differ according to the channel considered, so the "Log Scale" option allows for a better reading and analysis of the run.

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

Once this step has been completed, the analysis can proceed for the 5 channels considered FAM, HEX, Texas Red, CY5 and Quasar 705.

**To validate the assay**, the Ct values for the controls must be the following (Table 4). Outside of these values, the experiment cannot be validated.

**Table 4: Run validation**

| Positive Control |                   |
|------------------|-------------------|
| FAM              | Ct ≤ 30           |
| HEX              | Ct ≤ 30           |
| Texas Red        | Ct ≤ 30           |
| Cy5              | Ct ≤ 30           |
| Quasar 705       | Ct ≤ 30           |
| Negative Control |                   |
| FAM              | Ct not determined |
| HEX              | Ct not determined |
| Texas Red        | Ct not determined |
| Cy5              | Ct not determined |
| Quasar 705       | Ct not determined |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

## 12. Data analysis and interpretation

For clinical samples, the following results are possible:

\* Ct threshold for sample: + Positive → positive Ct < 45

For SARS-CoV-2 targets and mutations in samples:

Analyse the results as described below.

**SARS-CoV-2 (Texas Red), E484K mutation (FAM), E484Q mutation (HEX), L452R mutation (Quasar 705), Endogenous internal control (Cy5) : Ct < 45**

Results are analysed in the **FAM, HEX, Texas Red, Quasar 705** and **Cy5** channels (Table 5).

**Table 5: Detection of SARS-CoV-2 and mutations of interest**

| PCR signal |     |                  |            |     | Test validity                                                                          |
|------------|-----|------------------|------------|-----|----------------------------------------------------------------------------------------|
| FAM        | HEX | Texas Red        | Quasar 705 | Cy5 |                                                                                        |
| +          | -   | +/- <sup>1</sup> | -          | +/- | SARS-CoV-2 positive and presence of E484K mutation                                     |
| +          | -   | +/-              | +          | +/- | SARS-CoV-2 positive and presence of E484K and L452R mutations                          |
| -          | +   | +/-              | -          | +/- | SARS-CoV-2 positive and presence of E484Q mutation                                     |
| -          | +   | +/-              | +          | +/- | SARS-CoV-2 positive and presence of E484Q and L452R mutations                          |
| -          | -   | +/-              | +          | +/- | SARS-CoV-2 positive and presence of L452R mutation                                     |
| +          | +   | +/-              | -          | +/- | SARS-CoV-2 positive and and presumed positive for E484K or E484Q mutation <sup>2</sup> |
| -          | -   | +                | -          | +/- | SARS-CoV-2 positive without mutations E484K, E484Q, L452R                              |
| -          | -   | -                | -          | +   | SARS-CoV-2 negative                                                                    |
| -          | -   | -                | -          | -   | NI <sup>3</sup>                                                                        |

<sup>1</sup>+/- : if a signal is detected in one of the FAM/HEX/Quasar 705 channels but not in the Texas Red channel, the sample is considered to be SARS-CoV-2 positive because one or more mutations of interest are revealed.

<sup>2</sup>Additional confirmatory analysis can be performed in case it is necessary to differentiate between the E484K and E484Q mutation.

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

<sup>3</sup>NI : not interpretable because of RT-PCR inhibition or extraction problems: no conclusion can be given. It is then recommended to collect a new sample and/or repeat the extraction and/or dilute the sample 5 times.

## **Limitations on use and interpretation:**

- ❖ The EurobioPlex SARS-CoV-2 SNPs kit is used for first-line diagnostic purposes.
- ❖ All samples should be treated as potentially infected with SARS-CoV-2, and local biosafety regulations should be carefully followed.
- ❖ Interpretation of results should consider the possibility of false negatives and false positives.

False negatives may be due to:

- Inadequate collection of samples, or incorrect storage,
- Samples outside the viremia phase,
- Inappropriate extraction conditions or use of non-validated PCR instruments,
- Experimentation not respecting all elements of these instructions for use.

False positives may be due to:

- Contamination due to mishandling of high positive samples, the positive control, or the PCR amplification products,
  - Failure to follow the procedure described in these instructions for use, especially to avoid sources of contamination.
- ❖ All results should be interpreted by medical professionals in the context of the patient's clinical situation, history and symptoms.
  - ❖ This assay does not exclude the presence of pathogens other than SARS-CoV-2.
  - ❖ A negative result of this assay does not absolutely exclude a possible SARS-CoV-2 infection.

## 13. Performance analysis

### **Detection limitation/analytical sensitivity on CP:**

SARS-CoV-2 : 5 copies/μL

E484K mutation : 5 copies/μL

E484Q mutation : 5 copies/μL

L452R mutation : 10 copies/μL

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

## Analytical detection/sensitivity limitation on a synthetic SARS-CoV-2 RNA of Beta variant (Twist Bioscience):

SARS-CoV-2 : 1 copy/ $\mu$ L

E484K mutation : 12,5 copies/ $\mu$ L

E484Q mutation : N/A

L452R mutation : N/A

NA : not applicable.

## Analysis of two reference samples provided by the CNR des Virus des Infections Respiratoires (Hospices Civils de Lyon) from viral culture supernatant\* :

### - Beta variant

SARS-CoV-2:  $5.10^{e3}$  copies/ $\mu$ L

E484K mutation:  $5.10^{e3}$  copies/ $\mu$ L

### - Delta variant

ID sample: hCoV-19/France/NAQ-HCL021079453301/2021

SARS-CoV-2:  $4.10^{e4}$  copies/ $\mu$ L

L452R mutation:  $4.10^{e4}$  copies/ $\mu$ L

\* RNA extracted from VERO cell cultures.

The validation of the performance was performed on:

| 114 negative samples            |                                                                                                        |
|---------------------------------|--------------------------------------------------------------------------------------------------------|
| 94 SARS-CoV-2 negative samples  | Characterised with EurobioPlex EBX-041 CE-IVD kit                                                      |
| 20 positive respiratory samples | RSV-A/B, Parainfluenzae, Rhinovirus, Metapneumovirus, Bocavirus, CoV NL63, CoV OC43, CoV 229E, CoV KU1 |

| 128 positive samples                                                                                               |                                                                                                                          |
|--------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|
| 18 Alpha variant samples (lineage B.1.1.7)                                                                         | Characterised by sequencing                                                                                              |
| 22 Beta variant samples (lineage B.1.351)                                                                          | Characterised by sequencing                                                                                              |
| 3 samples with E484K mutation                                                                                      | Characterised by sequencing, EurobioPlex EBX-041 CE-IVD kit and Eurofins Genomics ViroBOAR Spike RT-PCR (SARS-CoV-2) kit |
| 63 SARS-CoV-2 positive samples not tested for the E484Q and L452R mutations, and not containing the E484K mutation | Characterised with EurobioPlex EBX-041 CE-IVD kit and Seegene Allplex SARS-CoV-2 Variants I Assay                        |
| 14 SARS-CoV-2 positive samples with L452R mutation <sup>1</sup>                                                    | Characterised with EurobioPlex EBX-041 CE-IVD kit and/or Seegene Allplex SARS-CoV-2 Variants I ou Variants II            |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

|                                                                |                                                                |
|----------------------------------------------------------------|----------------------------------------------------------------|
| 5 SARS-CoV-2 positive samples with E484Q mutation <sup>2</sup> | Characterised with EurobioPlex EBX-041 CE-IVD kit              |
| 1 SARS-CoV-2 positive sample with L452R and E484Q mutations    | Characterised by sequencing and EurobioPlex EBX-041 CE-IVD kit |
| 2 SARS-CoV-2 positive samples with L452Q mutation              | Characterised by sequencing and EurobioPlex EBX-041 CE-IVD kit |

<sup>1</sup>Samples being sequenced for confirmation of the L452R mutation, or already characterised by sequencing as Delta variant (lineage B.1.617.2), or lineage A.27, or lineage C.36

<sup>2</sup>Samples being sequenced for confirmation of the E484Q mutation.

For this study, sample extraction was performed on Maelstrom 9600 (TANBead) with the OptiPure Viral Auto Plate kit (ref: W665A46), and RT-PCR on CFX96 (Bio-Rad).

## SARS-CoV-2 general performance

|                                        |                            | EBX-046: SARS-CoV-2 detection channel (Texas Red) |                     |
|----------------------------------------|----------------------------|---------------------------------------------------|---------------------|
|                                        |                            | SARS-CoV-2 positive                               | SARS-CoV-2 negative |
| Pre-tested<br>EBX-041 or<br>sequencing | <b>SARS-CoV-2 positive</b> | 128                                               | 0                   |
|                                        | <b>SARS-CoV-2 negative</b> | 0                                                 | 114                 |

SARS-CoV-2 positive samples detected by the EBX-046 kit and pre-tested with the EBX-041 CE-IVD reference kit have Ct values between 9 and 37.

**General sensitivity:** > 99% (128/128)  
**General specificity:** > 99% (114/114)  
**General concordance:** > 99% (242/242)

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

Negative samples pre-tested with the EBX-041 CE-IVD kit:  
specificity for mutations of interest

| (-) negative             |       | SARS-CoV-2<br>(RdRp1/RdRp2 targets) | Mutations |       |       |
|--------------------------|-------|-------------------------------------|-----------|-------|-------|
|                          |       |                                     | E484K     | E484Q | L452R |
| SARS-CoV-2 negative      | N=114 | - (114/114)                         |           |       |       |
| Specificity to mutations |       | 100%                                |           |       |       |

Positive samples pre-tested with the EBX-041 CE-IVD kit or by sequencing:  
sensitivity for SARS-CoV-2 and specificity for mutations of interest

| N=128<br>(+) positive / (-) negative                        |      | SARS-CoV-2<br>(RdRp1/RdRp2 targets) | Mutations |       |       |
|-------------------------------------------------------------|------|-------------------------------------|-----------|-------|-------|
|                                                             |      |                                     | E484K     | E484Q | L452R |
| SARS-CoV-2 samples without E484K, E484Q, L452R mutations    | N=63 | + (63/63)                           | -         | -     | -     |
| Sensitivity to SARS-CoV-2                                   |      | 100%                                | N/A       |       |       |
| Specificity to mutations                                    |      | N/A                                 | 100%      |       |       |
| Alpha variant samples without E484K, E484Q, L452R mutations | N=18 | + (18/18)                           | -         | -     | -     |
| Sensitivity to SARS-CoV-2                                   |      | 100%                                | N/A       |       |       |
| Specificity to mutations                                    |      | N/A                                 | 100%      |       |       |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

|                                              |      |             |             |         |           |
|----------------------------------------------|------|-------------|-------------|---------|-----------|
| <b>Samples with E484K mutation</b>           | N=25 | + (25/25)   | + (25/25)   | -       | -         |
| Sensitivity to SARS-CoV-2                    |      | <b>100%</b> | N/A         |         |           |
| Specificity to mutations                     |      | N/A         | <b>100%</b> |         |           |
| <b>Samples with L452R mutation</b>           | N=14 | + (14/14)   | -           | -       | + (14/14) |
| Sensitivity to SARS-CoV-2                    |      | <b>100%</b> | N/A         |         |           |
| Specificity to mutations                     |      | N/A         | <b>100%</b> |         |           |
| <b>Samples with E484Q mutation</b>           | N=5  | + (5/5)     | -           | + (5/5) | -         |
| Sensitivity to SARS-CoV-2                    |      | <b>100%</b> | N/A         |         |           |
| Specificity to mutations                     |      | N/A         | <b>100%</b> |         |           |
| <b>Sample with L452R and E484Q mutations</b> | N=1  | + (1/1)     | -           | + (1/1) | + (1/1)   |
| Sensitivity to SARS-CoV-2                    |      | <b>100%</b> | N/A         |         |           |
| Specificity to mutations                     |      | N/A         | <b>100%</b> |         |           |
| <b>Samples with L452Q mutation</b>           | N=2  | + (2/2)     | -           | -       | -         |
| Sensitivity to SARS-CoV-2                    |      | <b>100%</b> | N/A         |         |           |
| Specificity to mutations                     |      | N/A         | <b>100%</b> |         |           |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

**484K sensitivity:** > 99% (25/25)  
**E484K specificity:** > 99% (217/217)  
**E484K concordance:** > 99% (242/242)

**E484Q sensitivity:** > 99% (6/6)  
**E484Q specificity:** > 99% (236/236)  
**E484Q concordance:** > 99% (242/242)

**L452R sensitivity:** > 99% (15/15)  
**L452R specificity:** > 99% (227/227)  
*No cross with L452Q mutation detected in the Quasar 705 channel (L452R)*  
**L452R concordance:** > 99% (242/242)

## Signal variability on FAM, HEX, Texas Red, Cy5 and Quasar 705 channels

- Intra-experimental variability:

| Mean CV % |       |       |            |                             |
|-----------|-------|-------|------------|-----------------------------|
| E484K     | E484Q | L452R | SARS-CoV-2 | Endogenous internal control |
| 1,70%     | 1,73% | 1,11% | 1,08%      | 2,03%                       |

CV: coefficient of variation

- Inter-batch variability:

| Mean CV % |       |       |            |                             |
|-----------|-------|-------|------------|-----------------------------|
| E484K     | E484Q | L452R | SARS-CoV-2 | Endogenous internal control |
| 2,12%     | 2,00% | 4,93% | 2,90%      | 3,83%                       |

CV: coefficient of variation

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

## Diagnostic specificity

This validation included:

- 44 samples negative for all known coronaviruses, including respiratory panel samples characterised for other viruses to test for cross-reactions.
- 7 positive samples for CoV NL63 (n=3) or CoV OC43 (n=2) or CoV 229E (n=1) or CoV KU1 (n=1).

| Number of samples | Positive for                       | Coronavirus | Status SARS-CoV-2 EBX-046 |
|-------------------|------------------------------------|-------------|---------------------------|
| 1                 | Adenovirus                         | Negative    | Negative                  |
| 2                 | RSV A ou B                         | Negative    | Negative                  |
| 7                 | Parainfluenzae                     | Negative    | Negative                  |
| 1                 | Rhinovirus                         | Negative    | Negative                  |
| 2                 | Metapneumovirus                    | Negative    | Negative                  |
| 1                 | Bocavirus                          | Negative    | Negative                  |
| 2                 | Legionella pneumoniae              | Negative    | Negative                  |
| 3                 | Mycoplasma pneumoniae              | Negative    | Negative                  |
| 3                 | Mycoplasma tuberculosis            | Negative    | Negative                  |
| 2                 | Bordetella pertussis               | Negative    | Negative                  |
| 1                 | Parechovirus Type 1                | Negative    | Negative                  |
| 1                 | Parechovirus Type 2                | Negative    | Negative                  |
| 1                 | Parechovirus Type 3                | Negative    | Negative                  |
| 1                 | Parechovirus Type 4                | Negative    | Negative                  |
| 1                 | Parechovirus Type 5                | Negative    | Negative                  |
| 1                 | Parechovirus Type 6                | Negative    | Negative                  |
| 3                 | Enterovirus                        | Negative    | Negative                  |
| 1                 | Streptococcus pyogenes             | Negative    | Negative                  |
| 3                 | Streptococcus pneumoniae           | Negative    | Negative                  |
| 2                 | Staph.epiderm/oralis/mitis/sanguis | Negative    | Negative                  |
| 3                 | Candida albicans                   | Negative    | Negative                  |
| 2                 | Pseudomonas aeruginosa             | Negative    | Negative                  |
| 3                 | CoV NL63                           | Positive    | Negative                  |
| 2                 | CoV OC43                           | Positive    | Negative                  |
| 1                 | CoV 229E                           | Positive    | Negative                  |
| 1                 | CoV KU1                            | Positive    | Negative                  |

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

In silico analysis was used and demonstrated the specificity of the primers and probes for SARS-CoV-2. For example, on the sequences of MERS and SARS viruses, no significant homology of SARS-CoV-2 primers and probes that could amplify these viruses was found.

## 14. Comparative analysis of the EBX-046 kit's diagnostic performance versus a commercial reference kit

The performance analysis of the EurobioPlex EBX-046 kit was conducted in comparison with another commercial CE-IVD reference kit by an external laboratory. The reference kit detects the presence of SARS-CoV-2 and targets the N501Y and E484K mutations, compared to our EBX-046 kit which detects the presence of SARS-CoV-2 and targets the E484K, E484Q and L452R mutations.

For this comparative study, 44 samples were tested and analysed. Four of the samples were retested with the EurobioPlex EBX-041 CE-IVD and Seegene Allplex SARS-CoV-2 Variants I Assay kits to validate their positivity to SARS-CoV-2.

|                              | EurobioPlex EBX-046 kit                |             |       |       | CE-IVD reference kit              |             |
|------------------------------|----------------------------------------|-------------|-------|-------|-----------------------------------|-------------|
|                              | SARS-CoV-2<br>(RdRp1/RdRp2<br>targets) | Mutations   |       |       | SARS-CoV-2<br>(RdRp/N<br>targets) | Mutations   |
|                              |                                        | E484K       | E484Q | L452R |                                   | E484K       |
|                              | + (42/44)                              | + (6/44)    | -     | -     | + (44/44)                         | + (6/44)    |
| Sensitivity to<br>SARS-CoV-2 | <b>95,45%*</b>                         | N/A         |       |       | <b>100%</b>                       | N/A         |
| Specificity to<br>mutations  | N/A                                    | <b>100%</b> |       |       | N/A                               | <b>100%</b> |

*\*Two samples were identified as negative for SARS-CoV-2. They have a very low viral charge, with no evident clinical significance (Ct>33,) and at the limit of sensitivity, which does not affect the possible use of the EBX-046 kit as a first-line test. This opinion, relating to weakly positive RT-PCR and the relevance of nasal swabbing, was issued by the Haut Conseil de la santé publique in November 2020 in an official press release (<https://www.hcsp.fr/Explore.cgi/AvisRapportsDomaine?clefr=963>).*

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## 15. Bibliography

Cherian, Sarah, Varsha Potdar, Santosh Jadhav, Pragya Yadav, Nivedita Gupta, Mousmi Das, Partha Rakshit, et al. 2021. « Convergent Evolution of SARS-CoV-2 Spike Mutations, L452R, E484Q and P681R, in the Second Wave of COVID-19 in Maharashtra, India ». Preprint. Molecular Biology. <https://doi.org/10.1101/2021.04.22.440932>.

Hu, Ben, Hua Guo, Peng Zhou, et Zheng-Li Shi. 2021. « Characteristics of SARS-CoV-2 and COVID-19 ». *Nature Reviews Microbiology* 19 (3): 141-54. <https://doi.org/10.1038/s41579-020-00459-7>.

Lauring, Adam S., et Emma B. Hodcroft. 2021. « Genetic Variants of SARS-CoV-2—What Do They Mean? » *JAMA* 325 (6): 529. <https://doi.org/10.1001/jama.2020.27124>.

Plante, Jessica A., Brooke M. Mitchell, Kenneth S. Plante, Kari Debbink, Scott C. Weaver, et Vineet D. Menachery. 2021. « The Variant Gambit: COVID-19's next Move ». *Cell Host & Microbe* 29 (4): 508-15. <https://doi.org/10.1016/j.chom.2021.02.020>.

Seyed Hosseini, Elahe, Narjes Riahi Kashani, Hossein Nikzad, Javid Azadbakht, Hassan Hassani Bafrani, et Hamed Haddad Kashani. 2020. « The Novel Coronavirus Disease-2019 (COVID-19): Mechanism of Action, Detection and Recent Therapeutic Strategies ». *Virology* 551 (décembre): 1-9. <https://doi.org/10.1016/j.virol.2020.08.011>.

Srivastava, Surabhi, Sofia Banu, Priya Singh, Divya Tej Sowpati, et Rakesh K. Mishra. 2021. « SARS-CoV-2 Genomics: An Indian Perspective on Sequencing Viral Variants ». *Journal of Biosciences* 46 (1): 22. <https://doi.org/10.1007/s12038-021-00145-7>.

European Centre for Disease Prevention and Control :  
<https://www.ecdc.europa.eu/en/covid-19/variants-concern>

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## 16. Quality control

In accordance with Eurobio's ISO EN 13485 certified quality management system, each batch of EurobioPlex SARS-CoV-2 SNPs is tested according to predefined specifications to ensure consistent product quality.

# Fiche technique EurobioPlex SARS-CoV-2 SNPs

---

## 17. Waste disposal

Eliminate all waste in accordance with the legislation on DASRI.

## 18. Incident report

Serious incidents related to the system must be reported to Eurobio Scientific and to the competent authority of the Member States in which the user and/or the patient is registered.

## 19. Technical assistance

For assistance with our products, please contact our technical support.

Eurobio Scientific's customer service can be contacted by e-mail at [adv@eurobio.fr](mailto:adv@eurobio.fr) or by phone at +33 (0)1.69.07.94.77.

7, avenue de Scandinavie  
91940 Les Ulis Cedex  
FRANCE