
eurobio

SCIENTIFIC

EurobioPlex

HDV qRT-PCR

REF EBX-004

24/48/96 réactions

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(p 2-19)

(p20-37)

(p38-55)

EurobioPlex

HDV qRT-PCR

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BIBLIOGRAPHIE

- (1) Rizzetto M, Canese MG, Arico S, Crivelli O, Trepo C, Bonino F, Verme G. Immunofluorescence detection of new antigen-antibody system (delta/anti-delta) associated to hepatitis B virus in liver and in serum of HBsAg carriers. *Gut* 1977; 18: 997-1003.
- (2) Brazas R, Ganem D. A cellular homolog of hepatitis delta antigen: implications for viral replication and evolution. *Science* 1996; 274:90-94.
- (3) Le Gal F, Gordien E, Affolabi D, Hanslik , Alloui C, De'ny P, Gault E. Quantification of Hepatitis Delta Virus RNA in Serum by Consensus Real-Time PCR Indicates Different Patterns of Virological Response to Interferon Therapy in Chronically Infected Patients. *J. Clin. Microbiol* 2005; 43: 2363–2369.
- (4) Zachou K, Yurdaydin C, Drebber U, Dalekos GN, Erhardt A, Cakaloglu Y, Degertekin H, Gurel S, Zeuzem S, Bozkaya H, Schlaphoff V, Dienes HP, Bock TC, Manns MP, Wedemeyer H; HDT-1 Study Group. Quantitative HBsAg and HDV-RNA levels in chronic delta hepatitis. *Liver Int.* 2010; 30:430-437.
- (5) EASL Clinical Practice Guidelines: Management of chronic hepatitis B virus infection. European Association for the Study of the Liver. *J. Hepat.* 2012; 57: 167–185
- (6) Le Gal F, Gault E, Ripault MP, Serpaggi J, Trinchet JC, Gordien E, Dény P. Eighth major clade for hepatitis delta virus. *Emerg Infect Dis.* 2006; 12: 1447-1450.
- (7) Chudy M, Development of the 1st International Standard for Hepatitis D Virus RNA –an update. 24th blood virology meeting, Ljubljana, 8-9 May 2013.
- (8) Manesis EK, Schina M, Le Gal F, Agelopoulou O, Papaioannou C, Kalligeros C, Arseniou V, Manolakopoulos S, Hadziyannis ES, Gault E, Koskinas J, Papatheodoridis G, Archimandritis AJ. Quantitative analysis of hepatitis D virus RNA and hepatitis B surface antigen serum levels in chronic delta hepatitis improves treatment monitoring. *Antivir Ther.* 2007;12: 381-388.
- (9) Fattovich G, Giustina G, Christensen E, Pantalena M, Zagni I, Realdi G, Schalm SW. Influence of hepatitis delta virus infection on morbidity and mortality in compensated cirrhosis type B. The European Concerted Action on Viral Hepatitis (Eurohep). *Gut.* 2000 Mar;46(3):420-6.
- (10) Radjef N, Gordien E, Ivaniushina V, Gault E, Anaïs P, Drugan T, Trinchet JC, Roulot D, Tamby M, Milinkovitch MC, Dény P. Molecular phylogenetic analyses indicate a wide and ancient radiation of African hepatitis delta virus, suggesting a deltavirus genus of at least seven major clades. *J Virol.* 2004 Mar;78(5):2537-44.
- (11) Castelnau C, Le Gal F, Ripault MP, Gordien E, Martinot-Peignoux M, Boyer N, Pham BN, Maylin S, Bedossa P, Dény P, Marcellin P, Gault E. Efficacy of peginterferon alpha-2b in chronic hepatitis delta: relevance of quantitative RT-PCR for follow-up. *Hepatology.* 2006 Sep;44(3):728-35.
- (12) Brichler S, Le Gal F, Butt A, Chevret S, Gordien E. Commercial real-time reverse transcriptase PCR assays can underestimate or fail to quantify hepatitis delta virus viremia. *Clin Gastroenterol Hepatol.* 2013 Jun;11(6):734-40
- (13) Brichler S, Le Gal F, Neri-Pinto F, Mansour W, Roulot D, Laperche S, Gordien E. Serological and molecular diagnosis of hepatitis delta virus infection: results of a French national quality control study. *J Clin Microbiol.* 2014 May;52(5):1694-7.

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
7, avenue de Scandinavie
91940 Les Ulis Cedex
France
+33 1 69 79 64 80

EurobioPlex *HDV qRT-PCR*

For **qualitative and quantitative** real time RT-PCR.

REF EBX-004

24/48/96 reactions

Version 2.03 – 25/11/2019

Validated only on CFX96™ real-time PCR detection system (Biorad) or Applied Biosystems® 7500 Real-Time PCR System (Applied Biosystems).

For any use of other equipment than CFX96 real-time PCR detection system or Applied Biosystems 7500 Real-Time PCR System, the results are not guaranteed

Storage conditions: keep all reagents between -15°C and -22°C until use and after first use **except for HDV standards and RC control which must be kept between -75°C an -85°C.**

Instructions for use

INTENDED USE

The HDV Eurobioplex test uses real time reverse transcription polymerase chain reaction (qRT-PCR) amplification and is designed for the qualitative and/or quantitative detection of hepatitis delta virus (HDV) from a nucleic acid RNA extract transcribed into cDNA by a reverse transcriptase (RT). The extract can come from plasma or serum of patients. This test is indicated to determine the viral load in patients during the natural history of the disease or treatment, as well as the possible replicative disposition of HDV infection. It especially targets patients infected by hepatitis B virus. The test is not intended for screening of banked blood or organs infected by HDV virus.

INTRODUCTION

Hepatitis D is an infection caused by the hepatitis Delta virus (HDV) which has been discovered at the end of the 70s by Mario Rizzetto. The HDV is a satellite virus of the hepatitis B virus (HBV). The HDV genome consists in a circular single-stranded RNA, pseudo-double strands, surrounded by two delta proteins (small p24 and large p27) forming ribonucleoprotein delta that is wrapped by the HBV envelope proteins or HBsAg. Consequently the HDV virus can only infect individuals being chronically infected with HBV.

HDV infection can manifest as a severe superinfection Delta often occurring in a patient already chronically infected by HBV, or as simultaneous acute co-infection by both viruses, the latter form being less severe. However the HDV is considered to be the virus responsible for the most serious liver disease, which is acute fulminating and chronic forms resulting in a faster progression to cirrhosis and hepatocellular carcinoma than in the mono HBV infection.

Approximately 10-15 million of individuals are infected worldwide by HDV. It is characterized by a very high genetic variability. Eight main genotypes (from 1 to 8) have been described with a characteristic geographical breakdown: the genotype 1 (HDV-1) is ubiquitous and is most frequently found; HDV-2 and -4 are found in Asia; HDV-3 in Northern parts of South America, and HDV-5 to-8 in sub-Saharan Africa. However due to the migration of populations these different genotypes can be found in many other countries.

Unlike HBV, HDV has no specific enzyme for replication of its genome, making an effective therapeutic strategy against this infection more difficult to establish. It is currently based only on the use of interferon alpha at high doses, and especially in its pegylated form. New anti- HDV specific drugs are currently being developed. It includes entry inhibitors and farnesylation inhibitors (key step of the viral cycle).

While the diagnosis of HDV infection is based on the detection of anti-delta antibodies, viral replication and monitoring of natural infection or antiviral treatment is based on the detection and quantification of the viral RNA by the technique of RT-PCR in real time in the blood of infected patients.

PRINCIPLE OF THE DETECTION

The Eurobioplex Hepatitis D test is a test of amplification of the ribonucleic acid (RNA) of the HDV virus as well as a control of RNA extraction and RT-PCR inhibition using one step RT-PCR amplification in real-time. Internal RNA control (RC) allows the control of all steps of the RNA extraction stage from blood samples (serum or plasma) through the stage of amplification by PCR in real time. It allows for checks of possible PCR inhibitors or technical problem.

The test is performed from the RNA extracted from the sample using a single reaction in a single well/tube. HDV virus is detected using a labeled probe FAM. Control of extraction and PCR inhibition is detected using a labeled probe HEX. All emit a specific fluorescence following its hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

This amplification system has been validated on a standard RNA containing a specific HDV virus insert (positive control HDV-QS), on the international standard from Paul Erlich Institute (1st WHO international standard for HD virus RNA for NAT testing WHO IS PEI: PEI code 7657/12), and on an internal RNA control (RC).

DESCRIPTION AND CONTENT OF THE KIT

The qPCR Hepatitis D kit is ready to use for the specific detection and quantification of the RNA of this virus. The kit contains reagents and enzymes for the amplification of Hepatitis D viral RNA and the internal control RNA (RC) by RT-PCR (Table 1). Fluorescence is emitted and recorded through optical measurements during the PCR using a fluorimeter by choosing the wavelength in the Table 2.

Table 1 : Note : well centrifuge the tubes before use

Component	Tubes -96 reactions	Tubes -48 reactions	Tubes -24 reactions	Description	Storage
RC	4 (white) tubes of 580 μ L	2 (white) tubes of 580 μ L	1 (white) tube of 580 μ L	Internal RNA Control (500 copies/ μ L) to add to the lysis buffer. For a manual/automated extraction: 10 μ L/sample for 70 μ L final volume of elution (Respect carefully the proportion depending on the volume of elution)	[-75°; -85°C] (defrost only once)

Component	Tubes -96 reactions	Tubes -48 reactions	Tubes -24 reactions	Description	Storage
Primers/Probes (Oligomix)	4 (orange) tubes of 35 μ L	2 (orange) tubes of 35 μ L	1 (orange) tube of 35 μ L	Primers and probes for <i>HDV</i> and RC	[-15°C ; -22°C]
HDV-QS (1)	5 (transparent) tubes of 1 standard of 20 μ L	3 (transparent) tubes of 1 standard of 20 μ L	2 (transparent) tubes of 1 standard of 20 μ L	Standard tube for quantification with: 1 quantity of <i>HDV</i> RNA synthetic matrix (1.10^7 copies/ 10 μ L) and 1 quantity of RC (100 copies/ μ L)	[-75°C ; -85°C] (defrost only once)
Molecular biology Water (CN-H2O)	4 (blue) tubes of 1.5 mL	2 (blue) tubes of 1.5 mL	1 (blue) tube of 1.5 mL	PCR grade water	[-15°C ; -22°C]
Reaction Mix (Enzyme)	4 (red) tubes of 370 μ L	2 (red) tubes of 370 μ L	1 (red) tube of 370 μ L	2x Reaction Mix containing Taq polymerase, 6 mM Mg-chloride and dNTPs	[-15°C ; -22°C]
Reverse transcriptase (RT)	2 (transparent) tubes of 13 μ L	1 (transparent) tube of 13 μ L	1 (transparent) tube of 7 μ L	Reverse transcriptase Enzyme	[-15°C ; -22°C]
	1	1	1	Instructions for use	

Table 2 :

Target	Fluorophore	Excitation	Emission
HDV	FAM	495 nm	520 nm
RNA Control (RC)	HEX	538 nm	555 nm

For any use with other equipment than CFX96™ real-time PCR detection system or Applied Biosystems® 7500 Real-Time PCR System, the results are not guaranteed

Alternatives:

- Channel **HEX** -> Channel VIC (Applied Biosystems 7500 Real-Time PCR System, Biorad CFX96™),

Required material not provided:

- ◇ Biological Hood
- ◇ qPCR instrument 96™ real-time PCR detection system or Applied Biosystems® 7500 Real-Time PCR System
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for qPCR
- ◇ Micropipettes
- ◇ Filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves

STORAGE

All reagents must be kept according to Table 1. Once prepared ARNs should be stored at -80°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing/defrosting cycles (> 3x for reagents stored between -15°C und -22°C) must be avoided, and could lead to decrease in sensitivity. Do not refreeze reagents at -75°C; -85°C.

Keep all reagents on ice during the experiment.

CAUTIONS AND NOTES

Read carefully instructions before starting.

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ Many freezing/defrosting cycles of the reagents must be avoided, and could lead to decrease in sensitivity.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ Define three working areas: 1) Isolation of RNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures. The use of gloves is mandatory. Appropriate methods of preparation of RNA for cDNA transcripts should be used.
- ◇ Always use filtered tips for micropipettes.
- ◇ Use specific lab coats and gloves in each working area.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

COLLECTION OF SAMPLES, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes containing EDTA anticoagulant for preparing plasma or without anticoagulant for serum. **Heparin is not usable due to its inhibitory effects on the PCR.**
- ◇ The samples must be immediately prepared and extracted or stored at +4°C up to 6h. To prepare the plasma or serum (once blood is coagulated), centrifuge blood 20 minutes at 800-1600 g and transfer the plasma or serum obtained in sterile tubes. If delays are longer only plasma and/or serum may be frozen from -20°C (1 week maximum) to -80°C (recommended if kept beyond one week).
- ◇ Transport of clinical samples must obey local regulations for this type of infectious agents. Human blood must always be considered as infectious and adequate protection procedures must be used.

Effects of biochemical parameters: very high concentrations of lipids (9 g/dL) or bilirubin (110 g/L) do not affect the results.

PROCEDURE

I- Extraction of RNA

RNA extraction kits are available from several manufacturers. You can use your own extraction homemade technique or a suitable commercial system by referring to the manufacturer's instructions. However, the equivalence between number of copies/microl and international units (IU/mL) have only been established for the following system:

m2000sp Abbott Automated extraction System Plasma vol: 500 µL Elution vol: 110 µL PCR sample vol: 10 µL

RC (read on the HEX channel) added upon extracting ensures that a negative result is not due to an extraction problem or due to the presence of too much inhibitors of RT-PCR. The RC solution concentration is 5000 copies/10 µL.

We recommend adding 10 µL of RC stock per sample to be extracted for a final elution volume of 70 µL. It is possible to add RC to lysis buffer for an automated extraction. Consequently, it will be distributed by the robot (the equivalent of 16 µL to 500 µL of sample per 110 µL elution volume for example with m2000sp automated extraction system, see above).

II- qRT-PCR

Prepare the reaction mix as follows with a surplus: number of test Nx120%.

	Vol (μL ; 1 test)
Reaction Mix (Enzyme)	12,5
Reverse Transcriptase (RT)	0,2
Primers/Probes (Oligomix)	1
Water (CN-H ₂ O)	1,3

To control that the PCR and extraction and real-time RT-PCR amplification steps worked properly, it is necessary to test:


For **qualitative analysis**:

- Hepatitis D positive control (CP): 10 μL of HDV-QS standard
- Negative control (CN): 10 μL of PCR water
- The control of extraction and amplification by RT-PCR (RC): 10 μL previously added to each sample before extraction (for a 70 μL elution)
- A positive control RC (CP-RC): 1.43 μL of RC + 8.57 μL of PCR water

Note: HDV-QS positive control contains a high concentration of matrix. Handling must be done carefully to avoid contamination.

For **quantitative analysis**:

In addition to the controls for qualitative analysis, it is necessary to prepare 4 additional dilutions of HDV-QS standard in addition to the tube supplied. The water supplied must be used as diluent. Assign the positive control HDV-QS supplied (10^7 copies/10 μL) as the highest standard in the first tube. Pipet 396 μL of sterile water in the tube 1, 36 μL in tube 2, 76 μL in tube 3, 16 μL in tube 4 and dilute as follows:

To generate a standard curve on a real-time PCR System, the 5 concentrations of standards must be used and defined as standard for HDV by specifying the corresponding number of copies/reaction or equivalent international unit (IU).


Example of plate layout for quantitative analysis

	1	2	3	4	5	6	7	8	9	10	11	12
A	HDV-QS1	Sample. 2										
B	HDV-QS2											
C	HDV-QS3											
D	HDV-QS4											
E	HDV-QS5											
F	CN											
G	CP-RC											
H	Sample. 2											

Sample: unknown RNA sample; HDV-QS: positive control/standard hepatitis D virus; CN: negative control (water); CP-RC: positive control RC

Protocole of qRT-PCR

- 1) Homogenize the Reaction Mix and centrifuge briefly.
- 2) Distribute the solutions according to the diagram below :

- 3) Close immediately with an adhesive film or caps to avoid contamination.
- 4) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 5) Program the qPCR instrument as follows:

Program	Temperature	Duration	Cycle(s)	
	46 °C	30 min	1	-
	95°C	2 min	1	-
Amplification	95°C	15 sec	45	-
	59°C	1 min		Fluorescence Acquisition

For any use with other equipment than CFX96™ real-time PCR detection system or Applied Biosystems® 7500 Real-Time PCR System, the results are not guaranteed

To generate a standard curve on a real-time PCR System, 5 HDV-QS standards and water (HDV=0) must be used and defined as a standard by specifying the corresponding number of molecules/reaction.

Input the concentrations of the HDV standards (FAM channel) during the setting of the thermocycler. The software will display a standard curve and will directly calculate the corresponding number of molecules by reaction. Depending on your quantity of DNA or biological fluid (ml) introduced (serum or plasma), this number of copies can be translated into a number of copies/μg of DNA or copies/mL of biological fluid.

The equivalence between the number of molecules HDV and the international standard has been validated for an RNA extraction systems (see table below). Thus, the number of IU/ml is directly posted on the thermocycler and the result in IU/mL of serum or plasma of patients is instantly calculated. There is a slight overestimation of the viral load for serum samples compared to plasma.

HDV-QS molecules/ 10 microlitres	m2000sp Abbott Automated extraction System Plasma vol: 500 μL Elution vol: 110 μL PCR sample vol: 10 μL
	Equivalence UI/ml
1E+07	3,40E+08
1E+05	3,40E+06
1E+04	3,40E+05
5E+02	1,70E+04
1E+02	3,40E+03

VALIDATION OF THE EXPERIMENT

Qualitative analysis:

- Negative control (CN): no amplification detected on the FAM channel. On the HEX channel, no amplification but it is possible to visualize a weak fluorescence signal. In this case the threshold bar must be positioned just above the maximum of the negative control signal.
- Positive control HDV-QS 1 (CP): PCR signal (FAM) at 15,8 Ct (+/- 3Ct)
- RC positive control: Ct on the HEX channel for the must be ≤ 38 .

Quantitative analysis:

- Negative control (CN): no amplification detected on the FAM channel. On the HEX channel, no amplification but it is possible to visualize a weak fluorescence signal. In this case the threshold bar must be positioned just above the maximum of the negative control signal.
- Ct values of HDV-QS standards on the FAM channel must be +/- 3 Ct of the values from the following table.

HDV-QS molecules/10 microlitres	C(t) FAM
1E+07	15,8
1E+05	21,8
1E+04	24,8
0,5E+03	29,8
1E+02	31,5

-- RC positive control: Ct on the HEX channel for the must be ≤ 38 .

- The efficiency of PCR on the FAM channel must be greater than 90% (equivalent to a "slope" $> -3,6$).

DATA ANALYSIS AND INTERPRETATION

Analysis of internal control RC in the samples :

If RC has been used according to the recommendations of the technical document, the expected value of Ct for the **RC should be ≤ 38** .

Two results can be obtained for RC:

1/ the RC control test is positive with $Ct \leq 38$: RNA has been properly extracted, and there is no inhibitor of RT-PCR. The result can be validated.

2/ the RC control test is negative or $Ct > 38$: either RNA was not properly extracted, or the RT did not function properly, or the presence of PCR inhibitors inhibits the PCR reaction. Then it is recommended to repeat the extraction or dilute the sample.

(Note: If the RC has not been used, no signal on the HEX channel is expected).

Qualitative analysis for HDV :

Positive (+): for Ct ≤ 45 (example curve a)

Negative (-): for Ct > 45 (example curve c)


Quantitative analysis for HDV :

Positive (+): for Ct ≤ 35 (example curve a)


Positive not quantifiable (PNQ): for $45 \geq Ct > 35 \Rightarrow$ indicates of a very low viral replication below the low threshold of quantification of the test (example curve b)

Negative (-): for Ct > 45 (example curve c)


Curve a : POSITIVE


Curve b : WEAKLY POSITIVE not quantifiable (PNQ)

Curve c : NEGATIVE

Or

In summary:

Possible results for clinical samples :

RC Ct valid : Ct≤38

RC Ct invalid : Ct>38

Depending on the value of RC Ct and the type of analysis (qualitative or quantitative), the HDV results will be usable or not

PCR Signal		Qualitative Analysis	Quantitative Analysis	Test validity/ comment
HDV (FAM)	RC (HEX)	Presence of HDV virus	Possible Quantification	
+ Ct≤35	Ct valid	Yes : positive sample	Yes	Valid for qualitative and quantitative analysis
+ 45≥Ct>35	Ct valid	Yes : positive sample	No : PNQ sample	Valid for qualitative but not quantitative analysis
+ Ct≤35	Ct invalid	Yes : positive sample	No	Valid for qualitative but not quantitative analysis * Possible inhibition of extraction or RT-PCR that prevents validation of the HDV quantification
+ 45≥Ct>35	Ct invalid	Yes : positive sample	No	Valid for qualitative but not quantitative analysis * Possible inhibition of extraction or RT-PCR that prevents validation of the HDV quantification
- Ct>45	Ct valid	No : negative sample	No : negative sample	Valid for qualitative analysis. No quantification since sample is negative
- Ct>45	Ct invalid	No possible interpretation	No possible interpretation	Not Valid for qualitative, nor quantitative analysis * Possible inhibition of extraction or RT-PCR

PNQ: weakly positive non-quantifiable

* If the expected information is not obtained, dilute 5x the sample. If the same result is obtained extract the sample again.

PERFORMANCE ANALYSIS


All experiments presented below were conducted on CFX96™ real-time PCR detection system (Biorad).

1- Analytical sensitivity and linearity for quantification:

The analytical sensitivity is defined as the smallest amount of target that can be precisely detected in a sample. This analytical detection limit was determined using HDV positive clinical samples in combination with a specific extraction method on the m2000sp automaton (Abbott Diagnostics with the RNA kit Abbott Diagnostic, m2000-RNA-Plasma-BA-500-70-V1 Protocol. Test sample of 500 microl and elution in 110 microl).

The detection limit is 1^{E+02} IU/ml. The number of International Units is derived from the HDV International Standard from Paul Erlich Institute (1st WHO international standard for HD virus RNA for NAT testing WHO IS PEI: PEI code 7657/12)


Sensitivity-Linearity of quantification: dilution series (1/5)

Linearity for quantification: validated from 1^{E+02} IU/ml to 1^{E+8,3} IU/ml for clinical sample.

2- Reproducibility within an experiment and between experiments :

N = 100 replicata for 2 concentrations of HDV international standard PEI ($1^{E+04,3}$ IU/ml and 1^{E+03} IU/ml) were carried out. The Coefficient of Variation within an experiment is 2.8% on log values for the 2 concentrations.

Reproducibility between experiments of the standard has been analyzed on the CFX96™ real-time PCR detection system (Biorad) (n = 110 measures): the median CV is 3%.

HDV-QS molecules/10 microlitres	Median C(t)	C(t) Std. Dev	Coefficient of variation inter runs (%) (n=13)
1^{E+07}	15,74	0,69	4,4
1^{E+05}	21,84	0,75	3,4
1^{E+04}	24,84	0,67	2,7
$0,5^{E+03}$	29,79	0,67	2,2
1^{E+02}	31,30	0,81	2,6

3-Specificity:

***Analytical specificity on non-HDV viremic samples:**

HDV negative plasma samples (negative samples for the detection of anti-HDV total antibodies) but viremic for other hepatotropic viruses were analyzed to determine the analytical specificity of the Kit Eurobioplex HDV, expressed as a negative result in the absence of the target. The analytical specificity was assessed on 105 specimens representing the following pathogenic viruses: HBV (n = 75), HIV (n = 10), HCV (n = 10), HAV (n = 5) or HEV (n=5).

No false positivity was found with the Eurobioplex HDV kits on the 105 samples tested.

4- Specificity and sensitivity :

***Sensitivity to the different HDV genotypes :**

8 major genotypes (HDV-1 to HDV-8) have been defined according to a nucleotide divergence within the genome of 20-35 %⁷⁻⁸. Each viral genotype emerged from a region and has now a defined geographic location. HDV-1 is ubiquitous, HDV-2 hits in Japan, Taiwan and in Yakutia (Russia); HDV-3 is exclusively observed in South America; HDV-4 in Taiwan and Japan, and HDV-5 to HDV-8 in sub-Saharan Africa. Owing to the migratory flows all these genotypes can potentially be found anywhere⁷⁻⁸. The performance of the HDV Eurobioplex has therefore been evaluated on 139 HDV positive RNA samples representing all viral genotypes and on 12 HDV negative RNA samples tested by the French reference national laboratory for HDV (CNR Delta-FNRL, Emmanuel Gordien, PhD.).

The following samples have been included in this part of the study: 33 HDV-1 Africa; 61 HDV-1 Europe/Asia; 1 HDV-2; 1 HDV-3; 1 HDV-4 (plasmid containing a copy of complete genome); 22 HDV-5; 7 HDV-6; 10 HDV-7 et 3 HDV-8. These 139 samples were quantified by the technique of CNR Delta and the kit Eurobioplex HDV EBX-004.

The 12 negative samples have been found negative using both techniques.

Eleven samples of very low viral loads (close to the thresholds of quantification) were detected by the 2 techniques (n = 3) or by one of the 2 techniques Eurobioplex HDV (n = 5) or CNR Delta (n = 3)

All other samples were detected and quantified equivalently by the two techniques. A coefficient of correlation $R^2=0,8$ was calculated indicating the excellent correlation of results between the 2 techniques whatever the viral genotype. The median value on all these samples was very close with 5.16 logUI/mL for the CNR Delta versus 5.05 logUI/mL for Eurobioplex HDV.

*** Global Sensitivity and specificity**

610 samples in total were used to assess the clinical/diagnostic sensitivity and specificity.

n=610		<i>Eurobio</i>	
		Pos	Neg
<i>CNR</i>	Pos	389	9*
	Neg	14*	198

* samples with very low viral load < 3logUI/mL

Clinical specificity : $198 / (198+14) = 93,4 \%$

Clinical sensitivity : $389 / (389+9) = 97,7 \%$
--

BIBLIOGRAPHY

- (1) Rizzetto M, Canese MG, Arico S, Crivelli O, Trepo C, Bonino F, Verme G. Immunofluorescence detection of new antigen-antibody system (delta/anti-delta) associated to hepatitis B virus in liver and in serum of HBsAg carriers. *Gut* 1977; 18: 997-1003.
- (2) Brazas R, Ganem D. A cellular homolog of hepatitis delta antigen: implications for viral replication and evolution. *Science* 1996; 274:90-94.
- (3) Le Gal F, Gordien E, Affolabi D, Hanslik , Alloui C, De´ny P, Gault E. Quantification of Hepatitis Delta Virus RNA in Serum by Consensus Real-Time PCR Indicates Different Patterns of Virological Response to Interferon Therapy in Chronically Infected Patients. *J. Clin. Microbiol* 2005; 43: 2363–2369.
- (4) Zachou K, Yurdaydin C, Drebber U, Dalekos GN, Erhardt A, Cakaloglu Y, Degertekin H, Gurel S, Zeuzem S, Bozkaya H, Schlaphoff V, Dienes HP, Bock TC, Manns MP, Wedemeyer H; HDT-1 Study Group. Quantitative HBsAg and HDV-RNA levels in chronic delta hepatitis. *Liver Int.* 2010; 30:430-437.
- (5) EASL Clinical Practice Guidelines: Management of chronic hepatitis B virus infection. European Association for the Study of the Liver. *J. Hepat.* 2012; 57: 167–185
- (6) Le Gal F, Gault E, Ripault MP, Serpaggi J, Trinchet JC, Gordien E, Dény P. Eighth major clade for hepatitis delta virus. *Emerg Infect Dis.* 2006; 12: 1447-1450.
- (7) Chudy M, Development of the 1st International Standard for Hepatitis D Virus RNA –an update. 24th blood virology meeting, Ljubljana, 8-9 May 2013.
- (8) Manesis EK, Schina M, Le Gal F, Agelopoulou O, Papaioannou C, Kalligeros C, Arseniou V, Manolakopoulos S, Hadziyannis ES, Gault E, Koskinas J, Papatheodoridis G, Archimandritis AJ. Quantitative analysis of hepatitis D virus RNA and hepatitis B surface antigen serum levels in chronic delta hepatitis improves treatment monitoring. *Antivir Ther.* 2007;12: 381-388.
- (9) Fattovich G, Giustina G, Christensen E, Pantalena M, Zagni I, Realdi G, Schalm SW. Influence of hepatitis delta virus infection on morbidity and mortality in compensated cirrhosis type B. The European Concerted Action on Viral Hepatitis (Eurohep). *Gut.* 2000 Mar;46(3):420-6.
- (10) Radjef N, Gordien E, Ivaniushina V, Gault E, Anaïs P, Drugan T, Trinchet JC, Roulot D, Tamby M, Milinkovitch MC, Dény P. Molecular phylogenetic analyses indicate a wide and ancient radiation of African hepatitis delta virus, suggesting a deltavirus genus of at least seven major clades. *J Virol.* 2004 Mar;78(5):2537-44.
- (11) Castelnau C, Le Gal F, Ripault MP, Gordien E, Martinot-Peignoux M, Boyer N, Pham BN, Maylin S, Bedossa P, Dény P, Marcellin P, Gault E. Efficacy of peginterferon alpha-2b in chronic hepatitis delta: relevance of quantitative RT-PCR for follow-up. *Hepatology.* 2006 Sep;44(3):728-35.
- (12) Brichler S, Le Gal F, Butt A, Chevret S, Gordien E. Commercial real-time reverse transcriptase PCR assays can underestimate or fail to quantify hepatitis delta virus viremia. *Clin Gastroenterol Hepatol.* 2013 Jun;11(6):734-40
- (13) Brichler S, Le Gal F, Neri-Pinto F, Mansour W, Roulot D, Laperche S, Gordien E. Serological and molecular diagnosis of hepatitis delta virus infection: results of a French national quality control study. *J Clin Microbiol.* 2014 May;52(5):1694-7.

WASTE ELIMINATION

Be in accordance with the law on the elimination of waste of clinical infectious material.

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
Keep protected from light

Manufacturer

CE marked product CE by the French notified body LNE GMED

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LITERATUR

- (1) Rizzetto M, Canese MG, Arico S, Crivelli O, Trepo C, Bonino F, Verme G. Immunofluorescence detection of new antigen-antibody system (delta/anti-delta) associated to hepatitis B virus in liver and in serum of HBsAg carriers. *Gut* 1977; 18: 997-1003.
- (2) Brazas R, Ganem D. A cellular homolog of hepatitis delta antigen: implications for viral replication and evolution. *Science* 1996; 274:90-94.
- (3) Le Gal F, Gordien E, Affolabi D, Hanslik , Alloui C, Dény P, Gault E. Quantification of Hepatitis Delta Virus RNA in Serum by Consensus Real-Time PCR Indicates Different Patterns of Virological Response to Interferon Therapy in Chronically Infected Patients. *J. Clin. Microbiol* 2005; 43: 2363–2369.
- (4) Zachou K, Yurdaydin C, Drebber U, Dalekos GN, Erhardt A, Cakaloglu Y, Degertekin H, Gurel S, Zeuzem S, Bozkaya H, Schlaphoff V, Dienes HP, Bock TC, Manns MP, Wedemeyer H; HDT-1 Study Group. Quantitative HBsAg and HDV-RNA levels in chronic delta hepatitis. *Liver Int.* 2010; 30:430-437.
- (5) EASL Clinical Practice Guidelines: Management of chronic hepatitis B virus infection. European Association for the Study of the Liver. *J. Hepat.* 2012; 57: 167–185
- (6) Le Gal F, Gault E, Ripault MP, Serpaggi J, Trinchet JC, Gordien E, Dény P. Eighth major clade for hepatitis delta virus. *Emerg Infect Dis.* 2006; 12: 1447-1450.
- (7) Chudy M, Development of the 1st International Standard for Hepatitis D Virus RNA –an update. 24th blood virology meeting, Ljubljana, 8-9 May 2013.
- (8) Manesis EK, Schina M, Le Gal F, Agelopoulou O, Papaioannou C, Kalligeros C, Arseniou V, Manolakopoulos S, Hadziyannis ES, Gault E, Koskinas J, Papatheodoridis G, Archimandritis AJ. Quantitative analysis of hepatitis D virus RNA and hepatitis B surface antigen serum levels in chronic delta hepatitis improves treatment monitoring. *Antivir Ther.* 2007;12: 381-388.
- (9) Fattovich G, Giustina G, Christensen E, Pantalena M, Zagni I, Realdi G, Schalm SW. Influence of hepatitis delta virus infection on morbidity and mortality in compensated cirrhosis type B. The European Concerted Action on Viral Hepatitis (Eurohep). *Gut.* 2000 Mar;46(3):420-6.
- (10) Radjef N, Gordien E, Ivaniushina V, Gault E, Anaïs P, Drugan T, Trinchet JC, Roulot D, Tamby M, Milinkovitch MC, Dény P. Molecular phylogenetic analyses indicate a wide and ancient radiation of African hepatitis delta virus, suggesting a deltavirus genus of at least seven major clades. *J Virol.* 2004 Mar;78(5):2537-44.
- (11) Castelnau C, Le Gal F, Ripault MP, Gordien E, Martinot-Peignoux M, Boyer N, Pham BN, Maylin S, Bedossa P, Dény P, Marcellin P, Gault E. Efficacy of peginterferon alpha-2b in chronic hepatitis delta: relevance of quantitative RT-PCR for follow-up. *Hepatology.* 2006 Sep;44(3):728-35.
- (12) Brichler S, Le Gal F, Butt A, Chevret S, Gordien E. Commercial real-time reverse transcriptase PCR assays can underestimate or fail to quantify hepatitis delta virus viremia. *Clin Gastroenterol Hepatol.* 2013 Jun;11(6):734-40
- (13) Brichler S, Le Gal F, Neri-Pinto F, Mansour W, Roulot D, Laperche S, Gordien E. Serological and molecular diagnosis of hepatitis delta virus infection: results of a French national quality control study. *J Clin Microbiol.* 2014 May;52(5):1694-7.

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