# **EurobioPlex**

## *HTLV1/HTLV2*

### **PCR EN TEMPS REEL**

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7, avenue de Scandinavie  
91953 Les Ulis Cedex  
FRANCE

**EurobioPlex**  
*HTLV1/HTLV2*  
**REAL-TIME PCR**

For **qualitative and quantitative** real-time PCR

**REF** EBX-008

**24/48/96 reactions**

**RUO**

Version 3.0 of 2019/04/24

**Validated on:**

- CFX96™ Real Time PCR detection system (Biorad) with analysis on CFX Manager v 3.1 (Biorad)
- LightCycler®480 Instrument II (Roche) with analysis on LightCycler® 480 software v1.5 (Roche)
- Applied Biosystems® 7500 Real-Time PCR Systems (Applied Biosystems) with analysis on 7500 Software v2

**Storage conditions:**

Keep all reagents between -15°C and -22°C until use and after first use

Instructions for use

## TABLE OF CONTENTS

<b>Intended Use .....</b>	<b>21</b>
<b>Introduction .....</b>	<b>21</b>
<b>Principle of detection .....</b>	<b>22</b>
<b>Description and content of the kit .....</b>	<b>23</b>
<b>Storage .....</b>	<b>24</b>
<b>Cautions and notes .....</b>	<b>24</b>
<b>Samples collection, transport and storage .....</b>	<b>25</b>
<b>Procedure .....</b>	<b>25</b>
<b>I- DNA extraction .....</b>	<b>25</b>
<b>II- Real-time PCR procedure .....</b>	<b>25</b>
II-1/ Qualitative PCR .....	26
II-2/ Quantitative PCR .....	26
II-3/ Protocol of real-time PCR .....	27
Diagram of the procedure .....	27
Detailed procedure .....	28
<b>Validation of the experiment .....</b>	<b>29</b>
<b>Data Analysis and Interpretation .....</b>	<b>30</b>
<b>Quantitative analysis.....</b>	<b>31</b>
<b>Performance analysis.....</b>	<b>32</b>
<b>Bibliography .....</b>	<b>34</b>
<b>Waste disposal .....</b>	<b>34</b>
<b>Symbols .....</b>	<b>35</b>

## INTENDED USE

The Eurobioplex HTLV1/HTLV2 test uses real-time polymerase chain reaction (PCR) amplification and is designed for the qualitative and quantitative detection of provirus HTLV1 and HTLV2 from a DNA nucleic acid extract or from a messenger RNA extract transcribed into cDNA. This test is indicated to detect and specifically differentiate the DNA of these proviruses, and thus to determine the occurrence of infection by HTLV1 and HTLV2 in humans, or to complete a proven or undetermined presence. Extracted DNA is the starting material for the Eurobioplex HTLV1-HTLV2 kit.

The Eurobioplex HTLV1-HTLV2 has been validated on the following specimen:

- Peripheral blood mononuclear cells (PBMC),
- cerebrospinal fluid,
- lymph node,
- skin,
- soft tissues.

## INTRODUCTION

HTLV1 is a retrovirus isolated in 1980 by the team of R. Gallo, from T cells of an adult patient with a Leukemia/Lymphoma T. This virus belongs to the family of the *Retroviridae*. Based on phylogenetic analysis of viral sequences, seven molecular subtypes of HTLV1 have been identified, in relation to specific endemic geographical areas. HTLV1 is essentially transmitted from cell to cell. Different cell types can be infected by the virus, especially CD4 T lymphocytes, but also CD8 T and B-lymphocytes, dendritic cells and macrophages. The multiplication of the virus is mainly linked to its clonal cell expansion, contrary to viruses such as HIV for which the expansion is due to non-productive viral cycles. Studies have estimated that 10 to 20 million people are infected with HTLV1 virus in the world. In Europe, HTLV1 infection is sporadic, corresponding to cases of importation (0.003%). In contrast, infection is endemic in the southwest of the Japan, South America (Brazil, Colombia, French Guiana, Surinam) and Central America, in the Caribbean basin (including the French West Indies), in Central and Eastern Africa. Thus, the prevalence of HTLV1 infection was evaluated between 1-5% among adults living in endemic areas: 1.2% in Togo, 1.76% in the Brazil, and 3% in the West Indies and French Guiana. The modes of contamination of HTLV1 virus in endemic areas are as follows:

- Sexual: mainly from men to women,
- Mother to child during long breast feeding (6 months),
- Blood transfusions (erythrocyte concentrates, platelet pellets) and intravenous drug abuse.

More often, HTLV1 infection remains asymptomatic. However, 3-10% of infected individuals develop severe disease, mainly:

- T Leukemia/lymphoma in the adult in 1-5% case. It corresponds to a monoclonal proliferation of activated CD4 T cells, and its modes of presentation are variable: leukemia (5%), lymphomatous (25%), indolent or chronic slow-onset (20%)
- Myelopathy: tropical spastic paraparesis / HTLV1 associated myelopathy (TSP/HAM) in 2-4% of cases

- Infectious dermatitis determined by a chronic exudative eczematous eruption,
- Intermediate uveitis of the adult; Sjögren syndrome
- Serious opportunistic infections such as anguillulosis, scabies and the dermatophyties.

However, HTLV1 infection may be responsible for other visceral complications. Other polymorphic manifestations related to infection by HTLV1 have been reported, such as lymphocytic alveolitis or diffuse centrilobular pulmonary nodules, dry eye or mouth syndrome, periodontitis, subcutaneous nodules of the hands and the olecranon, corticosteroid sensitive chronic arthritis, inflammatory muscle diseases associated with polymyositis/dermatomyositis. The prevalence of antibodies anti-HTLV1 is higher in patients with polymyositis/dermatomyositis than in the general population in Martinique (50% vs. 2.2%) and Japan (35% vs 11%).

The overall prevalence of neurological complications during infection by HTLV1 is evaluated between 2% and 4%. Myelopathy (TSP/HAM) is the neurological complication the best known and most common, resulting in motor disorders of the lower limbs characterized by a spastic paraparesis and genito-sphincter disorders (urinary frequency, mictional urgency).

Another HTLV human retrovirus is called HTLV2. It has been isolated for the first time in 1982 by R. Gallo from culture of splenic cells of a patient with malignant hemopathy considered initially as a variant of leukemia with tricholeucocytes (LT). HTLV2 endemic foci were observed in the Indian populations of South America in with a prevalence of 8 to 13 %. In Africa, cases of HTLV2 infection are present. In Europe, the prevalence in the general population is very low, but HTLV2 is present in France, as in other countries such as Italy and Spain, among intravenous drug users. In total, 4 genomic subtypes (A - D) are described to date. Infection with HTLV2 is mainly confined to CD8 cells and very weakly to CD4 T cells.

The detection of HTLV infection can be done by direct evidence of the virus either by viral culture or by viral genome screening using the technique of PCR. Anti-HTLV antibodies appear after a period of serological silence (serological window period) of on average 51 days after the counting. As with other retroviruses, carriers of antibodies are carriers of the virus, which persists throughout life. Antibodies are detected by ELISA and/or agglutination technique.

The detection by qPCR is sensitive, reproducible and brings an early detection of the presence of the virus.

## **PRINCIPLE OF DETECTION**

The Eurobioplex HTLV1/HTLV2 kit is a test using real-time amplification of viral DNA of these viruses (with differentiation), as well as a DNA extraction and PCR inhibition control. The test is performed from extracted DNA from a sample using a reaction in a single well/tube.

The DNA extraction and PCR inhibition control allows to check for variations that may occur during the DNA extraction step of biological samples and real-time PCR amplification. Thus, it ensures that a negative result may not be due to a bad DNA extraction, and/or the presence of PCR inhibitors in too large quantities.

This kit is composed of one triplex: HTLV1 + HTLV2 + control of DNA extraction and PCR inhibition (CI-PCR).

DNA of HTLV1 and HTLV2 is respectively detected using FAM and HEX labeled probes. DNA extraction and PCR inhibition control is detected using a CY5 labeled probe. All probes emit a specific fluorescence following their hydrolysis during the elongation of the amplification product. The measurement of the intensities of real-time fluorescence correlates with the accumulation of specific amplification products.

## DESCRIPTION AND CONTENT OF THE KIT

The real-time PCR Eurobioplex HTLV1-HTLV2 kit is ready to use for the specific detection and the quantification of these viruses.

Fluorescence is emitted and individually recorded through optical measurements during the PCR. The detection of the amplified fragment is performed by a fluorimeter using the channels shown in Table 1 below.

The kit contains reagents and enzyme for the amplification of HTLV1 and HTLV2 viral DNA, and the DNA extraction and PCR inhibition control (see Table 2).

**Table 1:**

Target	Fluorophore	Excitation	Emission
HTLV1	FAM	495 nm	515 nm
HTLV2	HEX	535 nm	555 nm
PCR inhibition control (CI-PCR)	Cy5	647 nm	667 nm

Equivalent channels on different real time PCR cyclers:

- Channel **FAM** (Systems ABI, SmartCycler II, Systems Mx, CFX96/Chromo4), Channel 510 (LC 480), Channel Green (RotorGene)
- Channel **HEX** (Chromo 4/CFX96, Mx Systems), Channel VIC (ABI Systems), Channel Alexa532 (SmartCycler II), Channel 580 (LC 480), Channel Yellow (RotorGene),
- Channel **Cy5** (CFX96/Chromo4, Systems ABI, Systems Mx), Channel Alexa647 (SmartCycler II), Channel 660 (LC 480), Channel Red (RotorGene)

Note: On LC480 instrument II, apply color compensation for the following channels: FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

**Table 2:**

Cap color	Components of the kit	96 reactions	48 reactions	24 reactions	Reconstitution
<b>Red</b>	ENZYME	4 x 375 µl	2 x 375 µl	375 µl	Ready to use
<b>Transparent</b>	OLIGOMIX	4 x 75 µl	2 x 75 µl	75 µl	Ready to use
<b>Yellow</b>	4 Standards: CP-1 : 3.10 <sup>6</sup> copies/µl CP-2 : 3.10 <sup>4</sup> copies/µl CP-3 : 3.10 <sup>3</sup> copies/µl CP-4 : 3.10 <sup>2</sup> copies/µl	4 x 75 µl ech CP -1 -2 -3 -4	2 x 75 µl ech CP -1 -2 -3 -4	75 µl ech CP -1 -2 -3 -4	Ready to use
<b>White</b>	PCR inhibition control (CI-PCR)	4 x 300 µl	2 x 300 µl	300 µl	Ready to use
<b>Blue</b>	Water = negative control (CN-H2O)	1 ml	1 ml	1 ml	Ready to use

Oligomix : contains the primers and probes of the triplex (HTLV1 + HTLV2 + CI-PCR)

**Required material not provided:**

- ◇ Biological Hood
- ◇ Real-time PCR instrument
- ◇ Micro centrifuge
- ◇ Vortex
- ◇ Plates / tubes for real-time PCR
- ◇ Micropipettes
- ◇ DNase-free RNase-free filter tips for micropipettes
- ◇ Sterile microtubes
- ◇ Gloves (powder free)

<b>STORAGE</b>
----------------

All reagents must be stored between -15 and -22°C.

All reagents can be used until the expiration date indicated on the kit.

Many freezing / defrosting cycles (> 3x) must be avoided, and could lead to decrease in sensitivity.

<b>CAUTIONS AND NOTES</b>
---------------------------

**Read carefully instructions before starting**

- ◇ The experiment must be performed by competent staff.
- ◇ Instruments must have been properly installed, calibrated and maintained according to the manufacturer's recommendations.
- ◇ Clinical samples are potentially infectious and must be processed under a laminar flow hood.
- ◇ The experiment must be performed according to good laboratory practices.
- ◇ Do not use this kit after expiration date.
- ◇ The kit is shipped with dry ice, and the components of the kits must arrive frozen. If one or more components are defrosted, or if the tubes have been damaged, contact Eurobio Scientific.
- ◇ Once defrosted, spin down briefly the tubes before use.
- ◇ It is recommended to define three working areas: 1) Isolation of DNA, 2) Preparation of the reaction mix and 3) Amplification / Detection of amplified products.
- ◇ Use specific lab coat and gloves (powder free) in each working area.
- ◇ Pipettes, reagents and other materials must not cross each area.
- ◇ Specific caution is required to preserve the purity of the reagents and reaction mixtures.
- ◇ Appropriate methods of preparation/extraction of DNA to produce high quality DNA and to be followed by a real-time PCR application should be used, particularly avoiding all sources of DNase contamination.
- ◇ Always use RNase-free DNase-free filtered tips for micropipettes.
- ◇ Do not pipette with mouth and do not eat, drink or smoke in the area.
- ◇ Avoid sprays.

## SAMPLES COLLECTION, TRANSPORT AND STORAGE

- ◇ Collect samples in sterile tubes. **The use of heparin as an anticoagulant is prohibited.**
- ◇ It is the responsibility of the user to master its own conditions of collection, transport and storage of samples, and extraction of DNA by suitable systems to produce DNA of good quality.
- ◇ The samples should be extracted immediately or stored following the recommendations in the table below (Table 3).

**Table 3 :**

Recommendations of maximum storage of samples before extraction	
Room temperature	2 h
4°C	16h
-20°C	Long term storage

- ◇ The user can refer to the World Health Organization or High Health Authority for storing samples.
- ◇ Transport of clinical samples must follow local regulations for this type of infectious agents.

## PROCEDURE

### **I- DNA extraction**

It is the user's responsibility that the extraction system used be compatible with downstream real time PCR technology. For this kit, we recommend to use extraction methods of viral DNA or of a messenger RNA extract transcribed into cDNA, from Peripheral blood mononuclear cells (PBMC), cerebrospinal fluid, lymph node, skin, or soft tissues, and refer to manufacturer's instructions of the extraction kit.

In the HTLV1-HTLV2 kit, CI-PCR on the CY5 channel can be added before extraction or in the PCR reaction. It ensures that a negative result is not due to an extraction problem or due to the presence of PCR inhibitors at high quantity.

We recommend the addition of 10 µl of CI-PCR per extraction. If the CI-PCR is added to control the PCR, CI-PCR is added to the reaction mix (1µl per PCR reaction). See real-time PCR protocol for details.

CI-PCR is available from Eurobio (Ref EurobioPlex EBX-002).

### **II- Real-time PCR Procedure**

#### **General comment:**

HTLV1-HTLV2 Standards and the DNA extraction and PCR inhibition control (CI-PCR) contain a very high concentration of matrix. Manipulations must be performed carefully to avoid contamination. To control the functioning of the PCR and steps of extraction and real-time PCR amplification, it is necessary to test a positive control, CP-1, as well as the negative control (water supplied = CN-H<sub>2</sub>O + CI-PCR) (see II-3/6) of real-time PCR procedure).

### **II-1/ Qualitative PCR**

No standard range is performed. Only positive control CP-1 and negative control are tested.


### **II-2/ Quantitative PCR**

It is necessary to perform a standard curve using the four ready to use standards tubes CP-1 to CP-4 in addition to a negative control. To generate a standard curve on a real-time PCR system, the concentrations of  $3.10^6$ ,  $3.10^4$ ,  $3.10^3$  and  $3.10^2$  copies/ $\mu$ l must be used and defined as standard for CP by specifying the corresponding number of copies/ $\mu$ l. In order to correct for possible experimental variation, we recommend testing CP standards in triplicates.

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### Diagram of the procedure

## Detailed Procedure

- 1) Homogenize the tube of Enzyme, and vortex Oligomix, CP standards and CI-PCR tubes before starting, and centrifuge.
- 2) Prepare a Mastermix as below; N is the number of PCR reactions. Plan to prepare enough reagents for at least N+3 reactions (refer to part 1-(a) or 1-(b) of the previous diagram according to the condition).

### Case (a): For samples extracted without CI-PCR

Number of reaction(s)	1	N+3
Enzyme	12.5 µl	(N+3) x 12.5 µl
Oligomix	2.5 µl	(N+3) x 2.5 µl
CI-PCR	1 µl	(N+3) x 1 µl
Total Volume Mastermix	16 µl*	(N+3) x 16 µl

### Case (b): For samples extracted with CI-PCR

Number of reaction(s)	1	N+3
Enzyme	12.5 µl	(N+3) x 12.5 µl
Oligomix	2.5 µl	(N+3) x 2.5 µl
Total Volume Mastermix	15 µl*	(N+3) x 15 µl

\* The volume difference between condition (a) or (b) has no effect on performance.

- 3) Homogenize mastermix prepared in 2) and centrifuge briefly.
- 4) Distribute 15 µl Mastermix\* with a micropipette and filtered tips in each tube or well of microplate for real time PCR.
- 5) Add 10 µl of DNA sample.
- 6) In parallel test the following controls:
  - Positive control(s):
    - o Qualitative test: 15 µl Mastermix + 10 µl CP-1
    - o Quantitative test: 15 µl Mastermix + 10 µl CP-1 to CP-4 to generate CP standard range
  - Negative control:
    - o Case (a): For samples extracted without CI-PCR
      - 15 µl Mastermix + 10 µl water supplied (CN-H2O)
    - o Case (b): For samples extracted with CI-PCR
      - 15 µl Mastermix + 9 µl water supplied (CN-H2O) + 1 µl CI-PCR
- 7) Close immediately with an adhesive film or transparent caps to avoid all contamination.
- 8) Centrifuge briefly to collect all the PCR reaction mix at the bottom of the tubes or plate.
- 9) Program the real-time PCR instrument as follows.

Program	Temperature	Duration	Cycle(s)	
<b>Denaturation</b>	95°C	3 min	1	-
<b>Amplification</b>	95°C	10 sec	40	-
	60°C	30 sec		Fluorescence Acquisition

Note 1: On LightCycler ® 480 systems (Roche), two optical systems are available: only "System II" is compliant with use of the kit. Apply color compensation for the following channels: FAM-HEX/VIC-Cy5 (465-510,533-580,618-660).

Note 2: For the Applied Biosystems systems, select "ROX" in "PASSIVE REFERENCE".


Note 3: On Rotorgene™, please calibrate the signal by clicking on "GAIN optimization".

Note 4: On CFX96 (Biorad), start the run using the v1.6 or later version of CFX Manager software, and analyse with v 3.1 (see § Validation of the Experiment)


## VALIDATION OF THE EXPERIMENT

The analysis of data post-acquisition on a CFX96 PCR instrument (Biorad) must be done with version 3.1 of CFX Manager Software (Biorad). In order to use this v3.1 version from a run started with an older version, follow the procedure below: at the end of the run, the data file with .pcrd suffix must be open and treated with version 3.1 of CFX Manager (Biorad).

If the run was done with CFX Manager v1.6 for instance, to open the data file with CFX Manager v3.1, click on CFX Manager v3.1 icon. The screen below appears.


- Click on File and select Open, then Data File

- Select the file you want to analyze and click on Open.

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Once this is done, the analysis can start.

**For the assay to be valid**, the Ct values for the controls must be the following (Table 4). Otherwise, the experiment is not valid.

**Table 4:**

Positive Control	
FAM	Ct ≤ 25
HEX	Ct ≤ 25
Negative Control	
FAM	Undetermined Ct
HEX	Undetermined Ct
CY5	Ct ≤ 35

Moreover, for quantitative PCR, in addition to the valid controls of Table 4, efficiency should be between 80 and 120% on the FAM and HEX channels, and the coefficient of correlation > 0.95.

## DATA ANALYSIS AND INTERPRETATION

### DNA extraction and PCR inhibition control in samples:

Two results can be obtained:

1/ the DNA extraction and PCR inhibition control (CI-PCR, channel Cy5) test is positive: DNA has been properly extracted, and there is no PCR inhibitor. The result can be validated.

2/ the DNA extraction and PCR inhibition control (CI-PCR, channel Cy5) test is negative: either DNA was not extracted, the PCR did not work well, or the presence of PCR inhibitors inhibits the PCR reaction. It is then recommended to repeat the extraction or dilute the sample, unless a specific signal appears in the targeted channels (FAM, HEX).

For clinical samples, the following results are possible:

\* Ct cut off for positivity of samples:

On all channels (FAM or HEX) : + Positive => Positive Ct (≤ 40)

PCR Signal			Presence of virus		Test validity/comment
FAM	HEX	CY5	HTLV1	HTLV2	
+*	-	+	Yes	No	<b>VALID</b>
-	+*	+	No	Yes	<b>VALID</b>
+*	+*	+	Yes	Yes	<b>VALID</b>
-	-	+	No	No	<b>VALID</b>
+*	+*	-	Yes	Yes	<b>VALID</b> Possible PCR inhibition or problem with extraction that does not prevent the detection of viruses.
+*	-	-	Yes	No possible interpretation	<b>VALID FOR HTLV1:</b> Possible PCR inhibition or problem with extraction that does not prevent the detection of HTLV1 <b>NOT VALID FOR HTLV2:</b> Possible PCR inhibition or problem with extraction that would prevent the detection of HTLV2 → <b>dilute first 5 x the sample and test again; if necessary redo an extraction</b>
-	+*	-	No possible interpretation	Yes	<b>VALID FOR HTLV2:</b> Possible PCR inhibition or problem with extraction that does not prevent the detection of HTLV2 <b>NOT VALID FOR HTLV1:</b> Possible PCR inhibition or problem with extraction that would prevent the detection of HTLV1 → <b>dilute first 5 x the sample and test again; if necessary redo an extraction</b>
-	-	-	No possible interpretation	No possible interpretation	Possible PCR inhibition or problem with extraction → <b>dilute first 5 x the sample and test again; if necessary redo an extraction</b>

## QUANTITATIVE ANALYSIS

Define the concentrations of the standards of HTLV1 (channel FAM) and HTLV2 (channel HEX) during the setting of the thermal cycler ( $3 \cdot 10^6$ ,  $3 \cdot 10^4$ ,  $3 \cdot 10^3$ , and  $3 \cdot 10^2$  copies/μl for CP-1, CP-2, CP-3 and CP-4, respectively). The software will draw the standard curve and will calculate directly the corresponding number of copies/μl for the samples.

Depending on the input quantity of DNA, volume of body liquid added, or number of peripheral blood mononuclear cells (PBMC), this number of copies/μl may be converted into a number of copy/μg (DNA) or copy/ml (biological fluid) or of copy/n cells (proviral load PVL/n cells).

It is also possible to normalize results by using a reference housekeeping gene of your choice by performing an independent PCR for this reference gene (not included in the kit). However, it is the responsibility of the user to validate such normalization system.

The results below quantification threshold cannot be quantitatively validated, but retain their quality of detectability (positivity validated by the fluorophore channel).

## PERFORMANCE ANALYSIS


### Limit of detection / Analytical sensitivity

100 copies/ $\mu\text{l}$  of positive control

### Linearity of quantification

From 300 copies/ $\mu\text{l}$  to  $3 \cdot 10^7$  copies/ $\mu\text{l}$  of positive control

### Example of experiment performed on real-time PCR thermocycleur CFX96 (Biorad):

### Validation on clinical samples

#### Study 1:

The tests were performed on CFX96 (Biorad).

The specificity and sensitivity were analyzed from DNA of 54 extracted samples from peripheral blood mononuclear cells (PBMC), cerebrospinal fluid, lymph nodes, skin, and soft tissues, from patients of different parts of the world and two independent reference centers.

➤ **Sensitivity:**

- HTLV1 positive: 31 PBMC samples, 7 cerebrospinal fluid samples, 1 lymph node sample, 1 skin sample and 2 other soft tissues,
- HTLV2 positive: 12 PBMC samples

**HTLV1: > 99%; HTLV2: 83%.**

➤ **Specificity:**

All of the tested negative samples are negative and concordant with pre-tests. **Specificity > 99% for HTLV1 and HTLV2.** No cross reaction between HTLV1 and HTLV2 was observed.

**Study 2:** This study has been performed from HTLV-1 positive whole blood samples. The tests were performed on ABI7500 (Applied Biosystems). The PBMC cells are separated from the blood samples (EDTA tubes) by centrifugation in ficoll density gradient. Then the QIAamp blood kit (Qiagen) was used.


➤ **Specificity: > 99%**

10 HIV positive samples and 10 zikavirus positive samples have been tested and found to be negative for HTLV-1.

➤ **Sensitivity: > 99%**

32 whole blood samples from subjects infected with HTLV-1, asymptomatic or suffering from paraparesis tropical spastic or myelopathy associated with HTLV-1 (TSP/HAM), have been tested, and all detected positive.

There was an excellent correlation between proviral loads obtained with Eurobioplex and the reference/pre-test method (graph below: yes = TSP/HAM; no = asymptomatic).

Legend: yes = TSP/HAM; no = asymptomatic

## BIBLIOGRAPHY

Lézin A, Amoureux MC, Fagour, Agbalika F and Césaire R, Quantitative real time PCR for HTLV1 and HTLV2, 17th International Conference on Human Retrovirology : HTLV and related viruses, 2015.

Dehée A, Césaire R, Désiré N, Lézin A, Bourdonné O, Béra O, Plumelle Y, Smadja D, Nicolas JC. Quantitation of HTLV1 proviral load by a Taqman real time PCR assay J Virol Methods. 2002 Apr;102 (1-2):37-51.

Gouhier E, Gaubert-Maréchal E, Abboud P, Couppié P, Nacher M. predictive factors of HTLV1-HIV coinfections in French Guiana. Am J Trop Med Hyg. 2013 Sep;89 (3):549-53.

de Mendoza C, Caballero E, Aguilera A, Pirón M, Ortiz de Lejarazu R, Rodríguez C, Cabezas T, González R, Treviño A, Soriano V; Spanish HIV-2/HTLV Group. HIV-2 and HTLV1 infections in Spain, a non-endemic region. AIDS Rev. 2014 Jul-Sep;16(3):152-9.

Paiva A, Casseb J. origin and prevalence of human T-lymphotropic virus type 1 (HTLV1) and type 2 (HTLV2) among indigenous populations in the americas. Rev Inst Med Trop Sao Paulo. 2015 Feb;57(1):1-14.

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## WASTE DISPOSAL

Be in accordance with the law on the elimination of waste of clinical infectious material.

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