

## Detection of anions in fire explosive materials

### Introduction:

In recent years, there have been cases of criminals using inorganic explosives such as ammonium nitrate, pyrotechnic agents, and black powder for criminal activities. By examining explosives and their explosive residues, the type of explosives and the nature of the case can be determined, providing direction and clues for investigation. At present, the main detection methods in China are chemical methods, which are cumbersome, low sensitivity, and have many interferences. Traditional methods have significant limitations. Ion chromatography, as a liquid phase branch, currently has multiple applications in environmental, food and other methods, and only requires one injection and can detect multiple anions and cations. The pre-treatment of the method is simple.

Detection items (Table 1):

Anion	Cl <sup>-</sup>	ClO <sub>3</sub> <sup>-</sup>	NO <sub>3</sub> <sup>-</sup>	SO <sub>4</sub> <sup>2-</sup>	ClO <sub>4</sub> <sup>-</sup>
-------	-----------------	-------------------------------	------------------------------	-------------------------------	-------------------------------

**Keywords:** Anions, Ion chromatography, Explosive.

### Instruments and equipment

- **Ion chromatograph:** CIC-D160<sup>+</sup>
- **Ultra pure water machine:** ECO-S15

Qingdao Shenghan Chromatograph Technology Co., Ltd



## Requirements

### Reagents

Unless otherwise specified, all reagents used are superior grade.  $\text{Cl}^-$ ,  $\text{ClO}_3^-$ ,  $\text{NO}_3^-$ ,  $\text{SO}_4^{2-}$  and  $\text{ClO}_4^-$  standard solution (1000 mg/L)

### Deionized Water

When preparing standard samples manually or diluting real samples, please use ASTM filtration and deionization requirements that meet the specifications listed in the table 2.

Table 2: Deionized water specification.

Specification	
Ions Resistivity	$\geq 18.25 \text{M}\Omega \cdot \text{cm}$
Organics-TOC	<10ppb
Iron/Transition Metals	<1ppb
Pyrogens	<0.03Eu/mL
Particulates (>0.2 $\mu\text{m}$ )	<1unit/mL
Colloids-Silica	<10ppb
Bacteria	<1cfu/mL

### Chromatography conditions:

Table 3: Analysis conditions

Instrument	CIC-D160 <sup>+</sup>
Eluent	7.0 mmol/L $\text{Na}_2\text{CO}_3$ +6.0 mmol/L $\text{NaHCO}_3$
Flow rate	0.5 mL/min
Injection volume	500 $\mu\text{L}$
Analytical Column	SH-AC-12C
Column oven temperature	35 $^\circ\text{C}$
Conductivity cell temperature	35 $^\circ\text{C}$
Suppressor current	35 mA

### Sample preparation

Extract the dust at the center of the explosion and extract a series of explosion dust in a certain distance, direction, and area from the explosion point. Simultaneously extract a series of blank control samples.

Weigh 1.5g of the sample into a beaker, dissolve it in a small amount of ultrapure water, transfer it to a 100mL volumetric flask, add ultrapure water to volume, shake well, take 1mL from it into a 500mL volumetric flask, add

ultrapure water to volume, shake well, After pass 0.22  $\mu\text{m}$  microporous membrane, analyze it on the machine.

### Standard chromatogram

Standard chromatogram, As shown in below:




Figure 1. Chromatogram of standard sample.

Table 4: Resolution of standard sample

No.	Anions	Concentration mg/L	Area/ $\mu\text{s}/\text{cm}^2\text{s}$	Resolution
1	$\text{Cl}^-$	1.0	582.34	6.22
2	$\text{ClO}_3^-$	1.0	187.02	3.73
3	$\text{NO}_3^-$	1.0	257.35	7.50
4	$\text{SO}_4^{2-}$	1.0	377.48	13.97
5	$\text{ClO}_4^-$	1.0	89.07	0.00

Table 5: Linearity of standard sample

Ions	Concentration/ mg/L	Linear equation	Linearity coefficient
$\text{Cl}^-$	0~1.0	$Y=5.544 \times 10^8 X + 1692000$	0.9991
$\text{ClO}_3^-$	0~1.0	$Y=1.866 \times 10^8 X - 1073000$	0.9998
$\text{NO}_3^-$	0~1.0	$Y=2.557 \times 10^8 X - 728200$	0.9996
$\text{SO}_4^{2-}$	0~1.0	$Y=3.658 \times 10^8 X + 1.128000$	0.9994
$\text{ClO}_4^-$	0~1.0	$Y=8.986 \times 10^7 X - 441200$	0.9997

### Minimum detectable concentration

Stable for 30 minutes after startup, baseline noise is 1.379ns/cm.

Calculate the theoretical minimum detection concentration with a signal-to-noise ratio of 3 times, and the theoretical minimum quantitative concentration with a signal-to-noise ratio of 10 times.

Table 6: Minimum detectable concentration

Ions	Theoretical	Theoretical
------	-------------	-------------

	minimum detection concentration	minimum quantitative concentration
Cl <sup>-</sup>	0.03777	0.1259
ClO <sub>3</sub> <sup>-</sup>	0.17997	0.5999
NO <sub>3</sub> <sup>-</sup>	0.16477	0.5492
SO <sub>4</sub> <sup>2-</sup>	0.13622	0.4541
ClO <sub>4</sub> <sup>-</sup>	0.19108	0.6369

### Sample chromatogram




Figure 2. Chromatogram of sample




Figure 3. Sample added scalar quantity 0.05mg/kg




Figure 4. Sample added scalar quantity 0.1mg/kg




Figure 5. Sample added scalar quantity 0.5mg/kg

### Results and calculations

Table 7: Sample test result

Ions	added scalar quantity mg/kg	Concentration mg/L	Recycling rates %
Cl <sup>-</sup>	0	0.1329	--
	0.05	0.1821	98.40
	0.1	0.2267	93.80
	0.5	0.6395	101.3
ClO <sub>3</sub> <sup>-</sup>	0	0.0829	--
	0.04	0.1238	102.3
	0.08	0.1593	95.50
	0.4	0.4724	97.38
NO <sub>3</sub> <sup>-</sup>	0	0.003009	--
	0.02	0.02221	96.01
	0.04	0.03931	90.75
	0.2	0.1895	93.25
SO <sub>4</sub> <sup>2-</sup>	0	0.03291	--
	0.02	0.04225	93.40
	0.04	0.05213	96.10
	0.2	0.1991	92.33
ClO <sub>4</sub> <sup>-</sup>	0	1.174	--
	0.5	1.687	102.6
	1	2.225	105.1
	2	3.445	113.6

Remarks: ① ND indicates not detected; ② The measured value has been deducted from the blank value;

③ There may be differences in testing results between different methods and laboratories;

### **Feasibility analysis and conclusion**

The above experiments prove that the detection method has good resolution and is suitable for the determination of the content of the components to be measured in the sample.

After testing, under chromatograph conditions, the five anions in explosives can be separated within 20 minutes, and the pre-treatment method is simple. Under these conditions, the stability of the tested ions is good, and the linear correlation coefficients can all be above 0.999, which meets the testing requirements and can be applied to the testing of fire and explosives in the public security industry.