

◆ Systems Ranging from Manual to Fully Automatic

- Fast steam distillation system
 - 3 - 5 minutes
- Automatic shutoff
- pH titration (optional)

◆ Applications

- Kjeldahl, TKN & Ammonia
- Alcohol in wine & beer
- Phenol & Fluoride



* Titrator sold separately

DigiPREP Steam Distillation Systems

FEATURES	DigiPREP DISTILLATION SYSTEM MODEL				
	100	200	300	400	500
CATALOG NUMBMER	010-035-001	010-035-002	010-035-003	010-035-004	010-035-005
Automatic Addition of H ₂ O		•	•	•	•
Automatic Addition of NaOH	•	•	•	•	•
Automatic Addition of H ₃ BO ₃				•	•
Programmable Reaction Time	•	•	•	•	•
Programmable Distillation Time	•	•	•	•	•
Automatic Steam Generation	•	•	•	•	•
Variable Steam Generating capacity (40% -100%)	•	•	•	•	•
Automatic Removal of Sample Waste by Suction			•	•	•
Number of Programs	1	1	10	99	99
Display Language Can Be Chosen by the User	•	•	•	•	•
Visual Error Reports	•	•	•	•	•
Acoustic Error Reports	•	•	•	•	•
Standby Operation Between Distillations	•	•	•	•	•
Level Monitoring for the Set of Canisters	•	•	•	•	•
Connection for Titration Unit					•
Voltage (volts)	230	230	230	230	230
Dimensions (W x H x D in mm)	410 x 675 x 410	410 x 675 x 410	410 x 675 x 410	410 x 675 x 410	410 x 675 x 410
Weight (kg)	34	34	35	35	36

CALL FOR MORE INFORMATION OR TO REQUEST A QUOTATION:

North America / International

Tel.: (800) 361-6820
+1 (514) 457-0701
Fax: (800) 253-5549
+1 (514) 457-4499

Europe

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

China

Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302

Brazil

Tel.: +55 (19) 8194-8488