


PCRBIO DNA Markers


- Ready to use
- Room temperature stable
- Quantitative



1.5 % TAE agarose gel

■ PCRBIOLadder I


PB40.11-01



1 % TAE agarose gel

■ PCRBIOLadder II


PB40.12-01



2 % TAE agarose gel

■ PCRBIOLadder III

PB40.13-01



1.7 % TAE agarose gel

■ PCRBIOLadder IV

PB40.14-01





PCRBIOSYSTEMS

simplifying research

Features

- Ready to use – load straight onto your gel
- Room temperature stable – store at 25°C
- Quantitative – helps to visualise PCR yield
- Wide range – 50bp to 10kb
- Evenly spaced bands
- Easy to identify reference bands

Environmental Note

Because the PCRBIOSYSTEMS Markers are room temperature stable we negate the need for cold temperature shipping. This reduces your carbon footprint to help maintain a healthy world, plus it saves you money.

Cat Number	Product Name	Pack Size
PB40.11-01	PCRBIOSYSTEMS Ladder I	100 lanes
PB40.11-05	PCRBIOSYSTEMS Ladder I	500 lanes
PB40.12-01	PCRBIOSYSTEMS Ladder II	100 lanes
PB40.12-05	PCRBIOSYSTEMS Ladder II	500 lanes
PB40.13-01	PCRBIOSYSTEMS Ladder III	100 lanes
PB40.13-05	PCRBIOSYSTEMS Ladder III	500 lanes
PB40.14-01	PCRBIOSYSTEMS Ladder IV	100 lanes
PB40.14-05	PCRBIOSYSTEMS Ladder IV	500 lanes

PCR Biosystems Ltd. | Aztec House | 397-405 Archway Road | London, N6 4ER

T: +44 (0) 203 930 8101 | E: info@pcrbio.com | Orders: sales@pcrbio.com

www.pcrbio.com