

Waste Water Set-Ups

Waste Water Set-Up AL400

Photometer AL400 with standard accessory, thermoreactor AL125, Infra-red data transmission module IRIM, tube stand, membrane filter set, instruction manual, warranty information

COD 0 - 150 mg/l and 0 - 1500 mg/l,
Ammonia 1 - 50 mg/l N,
Nitrate 1 - 30 mg/l N
Nitrite LR 0.01 - 0.3 mg/l N
Nitrogen 5 - 150 mg/l N
Phosphate 0.02 - 1 mg/l P / 0.06 - 3.5 mg/l PO₄

Waste Water Set-Up AL400

4214100



Waste Water Set-Up AL800

Spectrophotometer AL800, thermoreactor AL125, 5 round vials ø 24 mm, tube stand, membrane filter set, instruction manual, warranty information

COD 0 - 150 mg/l und 0 - 1500 mg/l,
Ammonia 1 - 50 mg/l N
Nitrate 1 - 30 mg/l N
Nitrite LR 0.01 - 0.3 mg/l N
Nitrogen 5 - 150 mg/l N
Phosphate 0.02 - 1 mg/l P / 0.06 - 3.5 mg/l PO₄

Waste Water Set-Up AL800

4712100



Reagents for waste water set-ups

COD 0-150 mg/l O ₂ (25 pc.), mercury free **	420710
(25 pc.)	420720
(150 pc.)	420725
COD 0-1500 mg/l O ₂ (25 pc.), mercury free **	420711
(150 pc.), mercury free **	420716
(25 pc.)	420721
(150 pc.)	420726
COD 0-15000 mg/l O ₂ (25 pc.), mercury free **	420712
(25 pc.)	420722
(150 pc.)	420727
** without chloride removal	
Ammonia VARIO HR tube test	4535650
Nitrate VARIO tube test	4535580
Nitrite LR VARIO powder pack	4530980
Nitrogen VARIO Total HR tube test	4535560
Phosphate VARIO Total HR tube test	4535210

Accessories for Waste Water Set-Ups

Set of round vials with lids Height 48 mm, Ø 24 mm	197629
Membrane filter set for use when preparing samples, 25 membrane filters 0.45 µm, 2 syringes 20 ml	366150
Vial stand for 6 round vials Ø 24 mm, acrylic glass	418951
Vial stand for 10 vials (Ø 16 mm or □ 13,5 mm), acrylic glass	418957
Automatic pipette*, 1 - 5 ml	419076
Pipette tips*, 1 - 5 ml (white), 100 pc.	419066
Automatic pipette**, 0.1 - 1 ml	419077
Pipette tips**, 0,1 - 1 ml (white), 1000 pc.	419073

* 0 - 150 mg/l and 0 - 1500mg/l ; ** 0 - 15000 mg/l